

INTELECT® RPW 2 GENERATION 2020

With 10 years of experience and 10,000 RPW devices sold all over the world Chattanooga® has learned a lot about what really matters to you and your patients and integrated it into our next generation clinical Radial Pressure Wave unit.

Smart Bluetooth module and the Intellect Connect App give the Intellect® RPW 2 device the ability to download data, produce reports and upload treatment data to additional devices within the practice. The Intellect Connect App will allow patient data to be collected and share between all Intellect Connect ready devices.

The new Intellect® RPW 2 combines streamlined design with ergonomics and patient-comfort with ease-of-use resulting into the most compact clinical unit currently on the market.

With "made in Germany & the USA" as a quality label we are proud to offer you the latest technology in a

FEATURES



SUGGESTED PROTOCOLS LIBRARY

Available at a finger click; treatment pictures, anatomical and pathological library, treatments guidelines.



LARGE CHOICE OF TRANSMITTERS*

17 different transmitters available for tendinopathies, trigger points, spine

OPTIONAL

V-ACTOR® HF*

Vibration frequency up to 50Hz for Vibration therapy and muscle relaxation.



LARGE GRAPHIC COLOR TOUCH SCREEN

High resolution 7" capacitive touch screen for greater customer experience.

WINDOWS 10 APP

The Intellect Connect App* gives great flexibility in the handling and processing patient data, allowing the export of patient treatment data and protocols from device to PC allowing not only the ability to produce reports but also to share treatments and patient data between multiple Intellect RPW 2 devices within a practice and in other smart enabled Chattanooga devices.



NEW FALCON HANDPIECE

allowing lower Bar pressures for an increased range of indications and applicators



EASY FIT FOR EVERY TREATMENT ROOM

Slim design 41 x 42 x 110cm (L x W x H), not a lot of space needed for this new RPW generation.

INTELECT® RPW 2 TRANSMIT-

A choice of standard transmitters for Radial Pressure Wave (RPW) indications. Each transmitter is unique by shape and material to deliver it's own energy flux density and wave profile.

TRIGGER POINTS / TENDINOPATHIES



Ro40

15 mm Energy beam transmitter with concave coupling surface, best for pain zones near the skin surface Penetration depth 0-35mm Intensity level: Medium



R15

15mm transmitter for any type of tendinopathies Penetration depth 0-35mm Intensity level: Medium



C15 CERAmA-x®

Ceramic Energy 15mm transmitter for any type of tendinopathies Penetration depth 0-35mm



A6

Acupuncture like transmitter 6mm for local trigger points Penetration depth 0-20mm Intensity level: Low



T10

Finger transmitter 10mm for deep target tendonopathies Penetration depth 0-40mm Intensity level: High



F15

White soft transmitter 15mm for superficial pain regions, muscles of mastication, cervical spine Penetration depth 0-30mm Intensity level: Very Low



ATLAS

Soft flexible tipped transmitter 15mm for highly sensitive regions and triggerpoints in the cervical spine Penetration depth 0-15mm Intensity level: Low

MYOFASCIAL TREATMENTS



D20-S

Standard Oscillator, 20mm transmitter for muscle and connective tissue Penetration depth 0-50mm Intensity level: Medium



D20-T

Golden Oscillator 20mm transmitter for professional athletes, High energy application Penetration depth 0-50mm Intensity level: Very High



D35-S

Large oscillator 35mm for large muscles and connective tissue Penetration depth 0-60mm Intensity level:

STANDARD ACCESSORIES

PART NUMBER (standard accessories)	DESCRIPTION
18080	FALCON® handpiece, including 2x projectile; 2x guiding tube; 1 cleaning brush, 1x wrench; 1x allen key
1x Power cord with EU plug for 2174KIT, 1x Power cord with North American plug for 2173KIT, 1x Power cord with Japanese plug for 2175KIT	
13-28660*	Quick start guide booklets (available in multiple languages)
22652	Transmitters tray
22651	Gel bottle support cup
4248	1x Transmission gel, 250 ml

SPINE TREATMENTS

Double Tip shape to allows treatments of the paravertebral region.



SPINE-ACTORS®

3 sizes span from 20 to 35mm

FASCIA TREATMENTS

Unique shape especially designed to dissolve clots in the Fascia structures.



PERI-ACTOR® KNUCKLE

For Medium-sized and large fasciae thoracolumbar fasciae, plantar fasciae, trigger bands, shoulder and neck region



PERI-ACTOR SCOOP

Treatment at or below the edges of bones : Front iliac crest, edge of the tibia



PERI-ACTOR® SCRAPER

Treatment of all fasciae trigger band



PERI-ACTOR SPHERE

Punctual treatment of fasciae: Plantar fasciae, Arch of foot, Hand, Trigger points, Cervical spine Fasciae

VIBRATION THERAPY



V10

10mm for superficia trigger points



V25

25mm for muscle massage



V40

40mm for large muscles groups

ORDERING DETAILS

CODE:	DESCRIPTION
2174KIT	Intellect RPW 2 230V unit set including standard accessories
2173KIT	Intellect RPW 2 120V unit set including standard accessories
2175KIT	Intellect RPW 2 100V unit set including standard accessories