

Electric Blood Collection /Donation/ Dialysis Chair

Features

Functions

- Functions, backrest & largest lifting
- CPR- quick release back part
- Back and leg can be adjustable, easy to choose chair position and bed position
- 560mm ergonomics seat height
- Upgraded armrest, can be right-left, front-back adjustable. Donor can control universal position

Accessories Details

- Optional genuine leather or supebic cover, with artisanal making, comfortable hand feeling.
- Ergonomics design removable pillow
- Armrest is R120° concave curve, avoid donor's arm sliding, ensure the donor's safety and comfort.
- Adjustable ABS tray table, with antiskid convex edge, tank for cup and cellphone's/5V optional. Very convenient for donors.
- Full SS IV pole, bearing 15kg

Important Quality Control

- * European 4 sector ergonomics design
- * Baosteel from Fortune Global 500
- * Environmental plastic can be sterilized in 100°C, tensile up to 30MPa
- * 11 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact
- * Panasonic Robotic ensure
- * 360° full smooth welding
- * American motor, with UL/ROHS/EN standard. Noise under 50dB, with lifetime 20,000times movement.



Description	Part No.
Blood Doner Chair-Electric	BD 2268
Size	L185XH(55-68)XW58cm

Electrical Blood Doner/Dialysis Chair

Functions

- 2 Functions, backrest & largest lifting
- CPR- quick release back part
- Back and leg can be adjustable, easy to choose chair position and bed position
- 560mm ergonomics seat height
- Upgraded armrest, can be right-left, front-back adjustable. Donor can control universal position

Accessories Details

- Optional genuine leather or supebic cover, with artisanal making, comfortable hand feeling.
- Ergonomics design removable pillow
- Armrest is R120° concave curve, avoid donor's arm sliding, ensure the donor's safety and comfort.
- Adjustable ABS tray table, with antiskid convex edge, tank for cup and cellphone's/5V optional. Very convenient for donors.
- Full SS IV pole, bearing 15kg

Important Quality Control

- * European 4 sections ergonomics design
- * Baosteel from Fortune Global 500
- * Environmental plastic can be sterilized in 100°C, tensile up to 30MPa
- * 11 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact
- * Panasonic Robotic ensure
- * 360° full smooth welding
- * American motor, with UL/ROHS/EN standard. Noise under 50dB, with life-time 20,000times movement.



Description	Part No.
Blood Doner Chair-Electric	BD 2271
Size	L185XH(55-68)XW58cm