

Shinhwa Medical is dedicated to contributing to the betterment of the health of mankind with love and frontier technology.



Spine-MT^{CORE} SHINHWAMEDICAL Co., Ltd.



Spinal Decompression & Mobilization
A Global Standard for Spinal Therapy



Indications

- Herniated nucleus pulposus
- Degenerative disc disease
- Scoliosis
- Joint problems
- Muscle imbalance

Effects

- Reduce pain
- Increase vertebral flexibility
- Reduce muscle spasm
- Relieve nerve compression
- Resorption of protruding discs

#20, Cheongwon-gil, Ilgwang-myeon, Gijang-gun, Busan Korea 46047
TEL : +82-51-724-6655 FAX : +82-51-723-6644 E-mail : 8614493@hanmail.net

Spinal Decompression & Mobilization
A Global Standard for Spinal Therapy


www.shinhwamedical.co.kr

The core of *spinal therapy*

Spine-MT^{CORE}



SHM SHINHWAMEDICAL Co., Ltd.




The core of **spinal therapy**

Spine-MT **CORE**

SHORT TREATMENT!
HIGH SUCCESS RATE!
CUSTOMIZED!

Core Therapy



CORE **01** Decompression

- Feedback System
- Comfortable, Effective



CORE **02** Correction

- Target Set-up
- Function : 1. Unilateral 2. Rotation 3. Level



CORE **03** Mobilization Exercise

- Intensity Set-up : 1. Low 2. Medium 3. High
- Angel Set-up : 5 Steps



CORE **04** Program

- Manual Mode : Sophisticated Treatments
- Special Mode : Various Automatic Treatments

Shinhwa Medical is dedicated to contributing to the improvement of the health of mankind with love and frontier technology.



Spine-MT **CORE** SHINHWHA MEDICAL Co., Ltd.



Cervical Target Angle

Explanation

- Treatment principle(flexion-distraction of Cox)
- Function to cure the cervical region to be treated through target angle adjustment
- Adjustable from 0 to 25°



Lumbar Target Angle

Explanation

- Treatment principle(flexion-distraction) of Cox(doctor of chiropractic)
- Function to cure the lumbar region to be treated through target angle adjustment
- Adjustable from 0 to + 25°



Cervical Unilateral application

Explanation

- Functions according to medical examination and treatment principle
- Functions to apply treatment according to the herniated direction of the disk
- Adjustable from 0 to 20°



Lumbar Unilateral application

Explanation

- Functions according to medical examination and treatment principle
- Functions to apply treatment according to the herniated direction of the disk
- Adjustable from 0 to 20°



Cervical Decompression therapy

Explanation

- Application of feedback system
- Application of exercise function
- Application of decompression function
- Application of automatic program



Lumbar Rotation application

Explanation

- Treatment principle of Cyriax (orthopedic surgeon)
- Treatment according to the type of escape of the disk
- Adjust the range of lumbar rotation
- Adjustable rotation stop position from 0 to 18 °



Thoracic & Lumbar Mobilization therapy

Explanation

- By mobilizing the upper body to the left and right in a mutually functional manner, the muscles and joints of the thoracic and lumbar spine are activated.



Lumbar Decompression therapy

Explanation

- Application of feedback system
- Application of exercise function
- Application of decompression function
- Application of automatic program