

Dialysis Chair

Blood Collection/ Phlebotomy Chair

Chair designed according to ergonomic human engineering Electric to achieve various positions with comfortable mattress, PU cover, waterproof Cushion: PU foam, Four casters and two with brake. Fixed Height , 2 Motors with Three Function.

Specification:

Dimension:190*930*600mm (L*W*H)

Backrest length:830mm,

Leg rest length:420mm.

Backrest adjustable from :-12°to 75°

Leg rest adjustable from :-80°to 12°

Armrest size:550*150*60mm(L*W*D)

Height of armrest to seat :190mm.

Seat Height: 580mm

Pillow size:400*230*80mm(L*W*T)

Caster diameter: 100mm

Safe maximum load: 200kg



Hand Rail



Double Armrest (Optional)



Hand Control

Ordering Details

CODE : BD 1500/2M

Description: Phlebotomy/Blood Doner Chair –Elec.



Trendelenburg Position(-12°)



Manual Dialysis/Phlebotomy/Blood Doner Chair

Main Features

Fixed Height Low carbon steel solid structure, appearance of high strength fiber material and multi-layer baking paint, with durability and easy to clean with comfortable mattress, PU cover, waterproof Cushion: PU foam Upholstery, Adjustable by gas spring, Four casters and two with brake,

Specification:

Dimension:190*930*600mm (L*W*H)

Backrest length:830mm,

Leg rest length:420mm.

Backrest adjustable from :-12°to 75°

Leg rest adjustable from :-80°to 12°

Armrest size:550*150*60mm(L*W*D)

Height of armrest to seat :190mm.

Seat Height: 580mm

Pillow size:400*230*80mm(L*W*T)

Caster diameter: 100mm



Ordering Details

CODE : BD 1500/M

Description: Phlebotomy/Blood Doner Chair