




- Ergonomic design
- Backrest and legrest synchronized adjustment by gas spring
- Equipped with two adjustable and removable PU armrests
- Electrostatic powder painted main structure
- Leather upholstery for easy cleaning and resistant to disinfections
- Fixed legs
- Lockable Castors (optional)



PU Adjustable Armrest



Gas Spring Lever

TECHNICAL DETAILS	
Width:	630 mm
Length:	1750 mm
Height:	1250 mm
Seat Surface Height:	52 mm
Safe Working Capacity:	180kg