

DIALYSIS CHAIR MANUAL

Dialysis chair is a specialized chair that is designed specifically for patients undergoing dialysis treatment. It is typically used in hospitals, clinics, and other healthcare settings, and is designed to provide a comfortable and safe seating option for patients during their dialysis treatment. The chair is adjustable, with features such as reclining backs, adjustable armrests, and footrests, to allow patients to find a comfortable position during treatment.

FEATURES

- I.V pole with headrest and foot step, 4" castor all with brakes
- The fabric is pressed PVC artificial leather the chair includes user friendly features, such as adjustment angle of electric ,CPR function, headrest, foot step, Trendelenberg position
- Optional: manual CPR, dinner table



ORDERING DETAILS

CODE:	DESCRIPTION
BD 1550M	Dialysis Chair Manual

DIALYSIS CHAIR ELECTRIC

Comfortable dialysis chair is ergonomically designed, which can better provide patients with more comfortable environment, reduce the psychological pressure of patients, and make the process of blood donation or dialysis more comfortable. Ideal for dialysis centers, blood donation centers and nursing homes.

SPECIFICATION

▪ Overall size	1950 x 900 x 580 - 780mm
▪ Back section length	830 ± 20mm
▪ Seat section length	520 ± 20mm
▪ Leg section length	530 ± 20mm
▪ Seat section width	600 ± 20mm
▪ Pillow size	380 x 250 x 70mm
▪ Armrest size	500 x 180 x 40mm
▪ Armrest swing out	0-180°
▪ Back section adjustment	-12° ± 75°
▪ Leg section adjustment	-75° ± 12°
▪ Height adjustment	580-780 ± 20mm
▪ Trendelenburg position	0-12°
▪ Loading capacity	240 kg

ORDERING DETAILS

CODE:	DESCRIPTION
BD 1550E	Electric Dialysis Chair

