

MEDICAL DENTAL CHAIR



Comfortable, contoured-round seat cushion. Surface material: PU
Adjustable torso support. Locks in place and easily height adjust for operators. Seat height adjustment. 450 -650 mm.
Lifting mode: steam-driven lift.
Chair foot: Zinc alloy + Polishing
Five legs cast steel base. Smooth running dual revolving castors.

ORDERING DETAILS

CODE:	DESCRIPTION
CHA 111	Chair with backrest and 1 Side Armrest

MEDICAL DENTAL CHAIR



Tilt adjustable, control/lock by lever switch underneath the seating
High density foam, Microfiber leather / PU leather /High Antimicrobial
Backrest: L360 x W300mm
Seat Size: 50cm
Seat: L450 x W410 x H80mm
Frame: Polished aluminum alloy

ORDERING DETAILS

CODE:	DESCRIPTION
CHA 104	Chair with backrest and 1 Side Armrest

PHYSICIAN CHAIR WITH BACKREST



Chrome arched base with heavy-duty, dual-wheel casters for smooth maneuverability. 360-degree swivel
Chrome heavy duty castors.
User Weight: 175kgs.
Size: 60 x 80 x 62cm

ORDERING DETAILS

CODE:	DESCRIPTION
CHA 103	Surgeon Stool with Backrest

THERAPIST CHAIR WITH BACKREST



Backrest: L340mm x W410mm
Seat: L420mm x W380mm,
Gas Lift: 140mm travel distance,
Base: Φ 560mm Aluminum alloy,
Castor: Φ 50mm, N.W.9 kg,
Upholstery: High density foam,
Cover: PVC / Customized

ORDERING DETAILS

CODE:	DESCRIPTION
CHA 101	Therapist Chair with Backrest