

CT707

Chiropractic Table

Decompression, Auto Drop and Auto Flexion IT Convergence Table



Product characteristic

Simplicity

Outstanding Design

Outstanding Structure

Electric mechanical Control

Precise manufacturing tech

Width	710mm
Length	1,845~2,245mm
Elevation range	590~890mm
Weight	135kg
Power	AC220V or 110V (Option) / 60Hz 150W
Permitted safety weight	120kg



Elevation



3 Auto Drop



Lumbar TracUon



Pelvis Auto---Flexion



Pelvis Angle



Cervial Traction



Manual Cervical Cox



Cervical Flexion Angle



Length control



CT 707

01M707 is IT Convergence spinal therapy eEquipment which has successfully evaluated by 30 professionals including doctors and combined all outstanding points of world famous chiropractic tables.



STEP-01
High quality
Auto 3Drop

Auto-cocking drop by electric motor and sensor(No compressor)

STEP-02
Cervical Traction
Fix the cervical
vertebrae with
the belt



STEP-03
Auto traction + Flexion
Circumduction

Memory Function by LCD monitor

STEP-04
Auto decompression
Extended step
by step



STEP-05
Manual cervical ext
Cervical extension

