

DERMA CHAIR WITH 4 MOTORS



High comfort bed equipped with 4 motors that control the height, the inclination of the backrest, the inclination of the seat and adjustment of the leg rest height movement. PU upholstery of high quality and easy cleaning. Standard hand remote is included to control all functions including the reset and lock features. The built-in foot pedal controls are located at the base on both sides and they control the height, backrest and footrest. The RESET button brings the chair to a seated position in the lowest setting to make it easier to get on and off the chair. The armrests can swing out of the way and be removed if desired. There is also a LOCK button that will disable all functions to prevent inadvertently moving the chair during a procedure.

SPECIFICATION

• Dimensions	: L185 x W83 x H60-88cm
• Weight Capacity	: 225kgs
• Headrest extension	: +/- 7.5cm
• Backrest angle	: 0° to 80°
• Seat angle	: 0° to 20°
• Footrest angle	: 0° to 85°
• Height up and down	: 60 to 88cm
• Footrest extension	: +/- 10.5cm
• Footprint	: 100 x 51cm
• Cushion thickness	: 10cm

FEATURES

- Built on that will disable all functions to prevent inadvertently moving the chair during a procedure. Specifications: x Built-in foot pedals allow control for 3 functions from both sides of the chair
- Tilt Function (Trendelenburg position capability)
- Width of the bed without armrest 62.5cm with armrest 83cm x Wide headrest that's extendable with a breather hole
- High quality vinyl upholstery for easy cleaning and maintenance
- Reset button to bring the chair to a seated position x Lock button to turn off all functions while performing a procedure
- Reversible and removable armrests

ORDERING DETAILS

CODE:	DESCRIPTION
DER 220D	Derma Chair with 4 Motors
DER 220B	Derma Chair with 3 Motors

DERMATOLOGY COUCH



High comfort bed equipped with 4 motors that control the height, the inclination of the backrest, the inclination of the seat and adjustment of the leg rest height movement. PU upholstery of high quality and easy cleaning. It's designed for small surgery, dermatology examination and plastic surgery. With the built-in foot pedals on both side, beautician can operate the chair hands free, it's not necessary to touch the hand remote which can meet the disinfection demand. (with safe lock function) full electric motion for backrest, angle, up down, and leg rest. 120 degree rotation system.

FEATURES

- Built-in foot pedals allow control for 3 functions from both sides of the chair
- Tilt Function (Trendelenburg position capability)
- Width : 62.5cm / Width with armrest 83 cm
- Wide headrest that's extendable with a breather hole
- High quality vinyl upholstery for easy cleaning and maintenance
- Reset button to bring the chair to a seated position
- Lock button to turn off all functions while performing a procedure
- Reversible and removable armrests

SPECIFICATION

▪ Dimensions	: L185 x W83 x H60-88cm
▪ Weight Capacity	: 225kgs
▪ Headrest extension	: +/- 7.5cm
▪ Backrest angle	: 0° to 80°
▪ Seat angle	: 0° to 20°
▪ Footrest angle	: 0° to 85°
▪ Height up and down	: 60 to 88cm
▪ Footrest extension	: +/- 10.5cm
▪ Footprint	: 100 x 51cm
▪ Cushion thickness	: 10cm

ORDERING DETAILS

CODE:	DESCRIPTION
DER 246EB/\$M	Dermatology Chair with 4 Motors



DERMATOLOGY COUCH




DER-19 is a reclining phlebotomy chairs and DERMA range of couches are Elegant and refined, to promote a classy image. Better still, they achieve practical positioning for almost any procedure. Perfect for plastic/cosmetic surgery, dermatology, maxifacial and ophthalmic procedures, the DERMA has two unique features: Separate headrest. Adjustable in height and tilt. 180° seat turn on both sides by single lever. Operating style" armrests, which clamp to the backrest rails. Adjustable in height angle and position on the rail. Suitable for arm/ breast surgery and aesthetics.

SPECIFICATION

- Safe working load: 1605kg
- Power height adjustment 63cm – 95cm
- Power backrest adjustment & power Trendelenburg tilt
- Manual assist leg rest extension
- Fully flat to upright positioning
- Separate adjustable headrest
- Hand remote control
- Total Length : 172cm
- 180° Seat adjusts manually
- Unit Weight : 80KG

ORDERING DETAILS

CODE:	DESCRIPTION
DER-19	Dermatology / Spa Couch



DERMATOLOGY COUCH



3 Section Electric Facial Treatment Chair is a versatile product that can save your business money, while also being extremely convenient to use. This treatment chair is split up into three different sections, with each being controlled by a localized motor. The motors control the treatment chair's height, backrest position and seat inclination. There is a foot rest also included with the treatment chair, which is controlled manually. This treatment chair can be set in the upright position for use as a chair or straightened out for use as a bed. Removable armrests are included. Another removable feature is the facial cushion which pops out to transform into a face rest.


STANDARD FEATURES

- Product Size L 180 x W 83 x H 57-83cm
- Standard Control - Handset Optional Foot Control
- Height - Adjustable, adjustable backrest, tilting, Trendelenburg, with leg rest, ergonomic
- Number of sections : 3
- Operation Mode: Electric
- Max Weight - 170kgs
- No. of Motors - 3
- G.W (lbs) - 80kg
- Voltage - 220v-240v
- Frequency - 50Hz-60Hz

ORDERING DETAILS

CODE:	DESCRIPTION
DER 224 A	Dermatology Chair with 3 Motor

ELECTRIC SPA COLUMN TABLE



Suitable for a range of treatments. The bed features a thick foam pad for superior comfort and three adjustable sections to make it easy to position clients for a variety of treatments. The bed has three built-in motors that control the height, backrest and central module movements. A simple handset is included to control the bed's movements. The design also includes double-arm adjustable supports, a removable headrest that covers a built-in facial hole and an extra postural cushion. There's even retractable wheels to make it easy to move the bed as needed. The bed is made to hold up to 385 pounds and is finished with high-quality PU leather upholstery for easy care and maintenance.

SPECIFICATION

- Product Size : L 202 x W 77 x H 56-85.8 cm
- Weight : (kg) 175kgs
- Motors : 3
- Voltage : 220V-2410V
- Frequency : 50Hz/60Hz Standard:
- Control : Handset Installation
- PU Leather

ORDERING DETAILS

CODE:	DESCRIPTION
DER 265	Spa Single Column Table with 3 Motor
DER 275	Spa Double Column Table with 3 Motor

DERMATOLOGY COUCH



DERMA range of couches are Elegant and refined, to promote a classy image. Better still, they achieve practical positioning for almost any procedure. Perfect for plastic/cosmetic surgery, dermatology, maxifacial and other procedures, the DER 221D is a multi-functional couch, ideal for spa body and face treatments, supplied with adjustable backrest and the headrest. It also features adjustable armrests that allow the customer to completely relax their neck and shoulder muscles in the Supine position.

SPECIFICATION

- Safe working load: 1605kg
- Power height adjustment 66 cm – 86cm
- Power backrest adjustment & power Trendelenburg tilt
- Manual assist hand rest function
- Fully flat to upright positioning
- Separate adjustable head rest
- Total Length : 194cm
- Width : 80cm including Armrests
- 3 Motors including 1 Column Actuators
- Hand Remote control
- Unit Weight : 80KG
- Voltage : 100-120V / 220-240V
- Frequency : 50Hz / 60Hz

ORDERING DETAILS

CODE:	DESCRIPTION
DER 221D	Dermatology / Spa Couch

DERMATOLOGY COUCH



DERMA range of couches are Elegant and refined, to promote a classy image. Better still, they achieve practical positioning for almost any procedure. Perfect for plastic/cosmetic surgery, dermatology, maxi facial and other procedures, The DER259 is a multi-functional couch, ideal for spa body and face treatments, supplied with adjustable backrest and face-hole headrest. It also features adjustable arm rests that allow the customer to completely relax their neck and shoulder muscles in the prone position.

SPECIFICATION

- Safe working load: 225kg
- Power height adjustment 63cm – 95cm
- Power backrest adjustment & power Trendelenburg tilt
- Manual assist hand rest function
- Fully flat to upright positioning
- Separate adjustable head rest
- Hand remote control
- Total Length : 198cm
- Width : 74cm
- 4 Motors including 2 Column Actuators
- Hand Remote control
- Unit Weight : 100KG
- Voltage : 100-120V / 220-240V
- Frequency : 50Hz / 60Hz

ORDERING DETAILS

CODE:	DESCRIPTION
DER259/4M	Dermatology / Spa Couch

DERMATOLOGY COUCH



The stylish facial treatment chair/massage bed is designed with a single column and four motors that control the height, backrest, footrest, and in clination of the entire bed. It features removable armrests, a headrest with a cutout for the face, and high-quality polyurethane padding that is easy to care for and clean.

FEATURES

- Perfect for massage therapists, medical spas, dermatology clinics, and plastic surgery practices
- Electric powered height, backrest, footrest, and Trendelenburg tilt adjustments
- Client max weight capacity 450 lbs
- Modern base design on a single column
- Armrests are reversible and removable
- Headrest adjusts to comfortably fit your patient and includes face hole to accommodate prone position
- Upholstery is extremely resistant to stains and regular wear and tear
- Pre-programmed memory function brings the chair back to upright position
- Includes hand remote
- Upholstery available in dark grey

ORDERING DETAILS

CODE:	DESCRIPTION
DER 222B	Dermatology Couch

DERMATOLOGY COUCH



DERMA range of couches are Elegant and refined, to promote a classy image. Better still, they achieve practical positioning for almost any procedure. Perfect for plastic/cosmetic surgery, dermatology, maxi facial and other procedures, The DER-D201 is a multi-functional couch, ideal for spa body and face treatments, supplied with adjustable backrest and the headrest. It also features adjustable armrests that allow the customer to completely relax their neck and shoulder muscles in the supine position.

SPECIFICATION

- Safe working load: 150 kg
- Power height adjustment 66cm – 86cm
- Power backrest adjustment
- Manual assist hand rest function
- Fully flat to upright positioning
- Separate adjustable headrest
- Total Length : 178 cm
- Width : 74 cm
- 3 Motors
- Hand Remote control
- Unit Weight : 65 KG
- Voltage : 100-120V / 220-240V
- Frequency : 50Hz / 60Hz

ORDERING DETAILS

CODE:	DESCRIPTION
DER-D201	Dermatology / Spa Couch