

Electric dynamic tilt table - TAB-120DM

Description:

- Electric dynamic tilt table mainly used Single motor control lift, single motor control bed stand.
- It use high-grade anti-pollution flame retardant PVC.
- Electric dynamic tilt table provides +85° tilting, adds to be used to diagnose dysautonomia or syncope.
- Electric dynamic tilt table can be used for long-term bedridden patients, to prevent the rise of orthostatic hypotension.
- This item also provide a Hi-Lo function for different situation,

Standard Specification:

- Electric dynamic tilt table have two electric motors, standard supports "Ti motion" brand for Hi-Lo & Tiling.
- Electric dynamic tilt table can lift for height variation and tilting by hand switch.
- Adjustable tiling featuring single motor, angle adjustable from 0° to +27°.
- Adjustable backrest featuring single gas struts, angle adjustable from 0° to +70°.
- Accessories: adjustable belts(3 pcs), adjustable foot plate, active table adjustable.
- Electric dynamic tilt table able to lift up to 175kg working weight from its lowest height.
- Retractable castors 100mm for excellent mobility and safety.
- Fire/stain/mildew/oil/water/wear resistant Med PVC available in 9 kinds of colors.

Dimensions:

- Upholstery Length: 249cm
- Upholstery Width: 75cm
- Adjustable Height: 67cm–103cm
- Foam Thickness: 5cm
- Foam Density: 36g/m2
Open: 249 x 75 x (67-103) cm
- Product Weight: 150/170kgs
- User Weight: 175kgs

