

CHIROPRACTIC TABLE



Ideal Plus electric chiropractic table has traction and can be used by people of different heights. Electrical control the bed height, good for chiropractors with different heights. This chiropractic massage traction table used for chiropractic wellness center, orthopaedics, Chinese medicine, chiropractic, sports injury, rehabilitation, sub-health studio.

Adjustable tension, Easy control levers, 4cm thick multi-layered sponge system, Adjustable angled headpiece.

Fire/stain/mildew/oil/water/wear resistant Med PVC available. This is Chiropractic drop table a versatile table providing flexibility for your practice. Density sponge, provide more comfortable environmental health. Lift for height variation by Foot switch or Hand switch. High quality chiropractic equipment and tables can offer, catering for all of your daily needs. Can be for the waist, neck, shoulders, foot and other parts of the role of adjuvant therapy to do the syndrome. Chiropractic drop table able to lift up to 220kg working weight from its lowest height. Our Chiropractic tables are high quality and easy to use.

SPECIFICATION

- Electric Height adjustment : 21.6"/55cm-38"/97cm
- Pneumatic Tilting
Headrest Angle Adjustment : -40° - 30°
- Pneumatic Leg rest
Angle Adjustment : -40° - 22°
- Leg plate length : 30.7"/78cm

OPTIONAL EXTRAS

Bandage, Adjustable footrest, Moving wheel transportation system, Stainless steel foot switch, Hand Controller, Paper Roll Holder

ORDERING DETAILS

CODE:	DESCRIPTION
FR3300	Chiropractic Electric Table

