



GYNACOLOGY CHAIR- WITHOUT TRAY



Standard configuration

GYN-246EB-1



Option : Guide rail & Foot salver version

GYN-246EB-2



Option : Guide rail version

GYN-246EB-3

• Standard configuration

- Built-in Foot Controller
- Reset
- Lock
- 175kg 386lbs MAX Max load bearing

• Options

- Hand Controller
- Heating
- Memory System
- Foot Controller
- Silicone Sponge/ Memory Foam

• Specifications

Max load bearing : 175 kg / 386 lbs Product size(L*W*H) : 1770x680x500 / 800 (mm)
Net weight :103 kg / 227 lbs Gross weight :116 kg / 256 lbs
Package size : 1540x810x605 (mm)≈0.75m³