

Medical Stool and Chair

The Medical stool for assistants comes with an adjustable footing and a rotating support arm. The seat height is adjustable. The padded seat is round and with the scuff resistant PU leather. These dental assistant stools are not only engineered to provide exceptional comfort, they are also ergonomically release musculoskeletal pain common in dentistry. Adjustable armrest support can easily move and locks, and can be converted for left or right-handed operators.

Key Features

- Ergonomics, comfort, durability
- Round seat cushion
- Full metal frame
- Highly stable 5-wheel aluminum alloy base

Specification

- Armrest size: L45 × W10 cm
- Seat size: $\Phi 40$ cm
- Gas lift cylinder: 160 (optional 180)
- Seat height: 45-65cm
- Variable height by gas lift cylinder
- Footing: $\Phi 43$ cm aluminum alloy
- Base: $\Phi 56$ cm aluminum alloy
- Castors: $\Phi 50$ mm silent castor
- Upholstery: High density PU foam
- Cover: PU leather / Microfiber leather



Ordering Details

CODE : CHA 112/BS

DESCRIPTION: Hight Adjustable Stool with Back Support



Microfiber Leather

				
MF6-Blue	MF7-Sky	MF8-Cyan	MF9-Green	MF10-Light Green
				
MF11-Brown	MF12-Curry	MF13-Light Brown	MF14-Wheat	MF15-Creamy White
				
MF16-Black	MF17-Charcoal Grey	MF18-Grey	MF19-Grey White	MF20-White

Ordering Details

CODE : CHA 112

DESCRIPTION: Hight Adjustable Stool