


Medical Stool and Chair

Specification

- Material: PU / Micro fiber
- The cushion size: 380mm
- Back of a chair size: 500*260mm
- Steam lever Adjust th: 120mm-180mm
- Lift way: Steam lifting
- Chair foot: Zinc alloy die casting
- Diameter of foot 540mm
- Load capacity : 550kg
- Caster: High quality Noiseless
- Seat height: 620mm
- Seat size: 480mm
- Packing weight: 13KG
- Packing size: 560*290*570mm



Ordering Details

CODE : CHA 018

DESCRIPTION: Hight Adjustable Stool with Back Support

CODE : CHA 018 /FR

DESCRIPTION: Hight Adjustable Stool with Back Support & Foot ring

Doctors Chair with Armrest & Backrest

Specifications

- Chair material: Microfiber leather, built-in high rebound sponge, plastic plate
- The cushion size: 445mm*400mm
- Back of a chair size: 400mm*260mm
- Armrest size: 270mm*95mm*45mm
- Chassis: Chassis back regulation 10
- Steam lever: Adjust the stroke 120 mm-180 mm
- Lift way: Steam lifting
- Chair foot material: Zinc alloy die casting
- Surface treatment: polishing
- Seat diameter: 560mm
- Load capacity: 550KG
- Castor: High quality Noiseless casters
- Armrest material: Super fibrous skin, built-in sponges injection plates
- Hghest seat cushion: 630mm
- Lowest seat cushion: 500mm
- Top of the armrest: 825mm
- Minimum support: 745mm
- Chair Weight: 15Kg
- Packing weight 16.5kg
- Packing size 560x570x290mm



Ordering Details

CODE : CHA 037

DESCRIPTION: Doctor Chair with Arm & Backrest

CODE : CHA 037A

DESCRIPTION: Doctor Chair with Arm & Backrest