


MAT TABLE GRAN COMFORT



Treatment table with electrically adjustable height control for Bobath/Mat and Vojta therapy.

FEATURES

- Electric motor with a thrust force of 8000 N
- Tubular protection for upholstery
- Class IM1 fire retardant padding and covering
- Transport wheels



SPECIFICATION

- Height adjustment range : 46 Δ 94 cm
- Safe working load : 300 kg
- Dimensions : 200 x 200 x 46 Δ 94 cm

OPTIONAL EXTRAS

Peripheral Control Bar 240029dp

ORDERING DETAILS

CODE:	DESCRIPTION
MAT 240029D	Mat Table Gran Comfort


MAT BOBATH PRO POWER



Treatment table with electrically adjustable height control for Bobath or Mat and Vojta therapy.

FEATURES

- Electric motor with a thrust force of 8000 N
- Tubular protection for upholstery
- Class IM1 fire retardant padding and covering
- Transport wheels



SPECIFICATION

- Height adjustment range : 46 Δ 94 cm
- Safe working load : 300 kg
- Dimensions : 200 x 120 x 46 Δ 94 cm

OPTIONAL EXTRAS

Peripheral Control Bar 240029dp

ORDERING DETAILS

CODE:	DESCRIPTION
MAT 240041D	Mat Bobath Pro Power

FIXED HEIGHT MAT / BOBATH TABLE - MAT 240032



Fixed Height Bobath table with steel frame for medical examination, therapy and massages.

FEATURES

- Bobath / Mat therapy treatment table
- Fire-retardant padding and covering class 1IM, washable

SPECIFICATION

- Safe working load : 160 kg
- Dimensions : 200 x 200 x 50 cm

OPTIONAL EXTRAS

Peripheral Control Bar 240029dp

ORDERING DETAILS

CODE:	DESCRIPTION
MAT 240032	Fixed Height Mat

PORTABLE ALUMINIUM MASSAGE COUCH



Portable Aluminum frame with powder coated paint. Adjustable back rest and light weight aluminum frame. Auto Cabling Lock System, Weight 14kg, user Weight 180Kg. Size : L188cm x W68 x H61 to 77cm.

ORDERING DETAILS

CODE:

POR AT 204

DESCRIPTION

Portable Aluminum Massage Couch

PORTABLE ALUMINIUM MASSAGE COUCH



Portable Aluminum frame with powder coated paint, Adjustable back rest and height, light weight aluminum frame. Auto Cabling Lock System, Weight 13kg, user Weight 180Kg. Size : L185cm x W70 x H61 to 82cm.

ORDERING DETAILS

CODE:

POR 263B

DESCRIPTION

Portable Aluminum Massage Couch

PORTABLE ALUMINIUM MASSAGE COUCH



Portable Aluminum frame with powder coated paint, Adjustable back rest and height. light weight aluminum frame. Auto Cabling Lock System, Weight 13kg, user Weight 180Kg. Size : L188cm x W68 x H61 to 77cm.

ORDERING DETAILS

CODE:

POR AT304

DESCRIPTION

Portable Aluminum Massage Couch

PORTABLE WOODEN MASSAGE COUCH



Portable Wooden frame with Adjustable back rest and height. Light weight Beach wood frame. 4cm thickness foam, Auto Cabling Lock System, Weight 15kg, user Weight 200Kg. Size : L188cm x W71 x H61 to 86cm.

ORDERING DETAILS

CODE:

POR WT304

DESCRIPTION

Portable Wooden Massage Couch

PORTABLE WOODEN MASSAGE COUCH



Portable Wooden frame with Adjustable back rest and height. light weight Beach wood frame. 5cm thickness foam, Auto Cabling Lock System, Weight 15kg, user Weight 200Kg. Size : L185cm x W70 x H62 to 87cm.

ORDERING DETAILS

CODE:

POR WT221B

DESCRIPTION

Portable Wooden Massage Couch

FIXED HEIGHT STAINLESS STEEL WITH PAINTED COUCH



Structure made by Steel Square Tube Frame. Epoxy Powder Coated, 5cm Thick Molded Foam Pad, Fire retardant upholstery, Adjustable Back Rest Assisted by gas spring & Paper Roll Holder. Size : L195cm x W 60cm X H 77cm. Weight Capacity : 200kgs.

ORDERING DETAILS

CODE:	DESCRIPTION
TAB 200 FHS	2 Section Fixed Height Couch

FIXED HEIGHT WOODEN COUCH



2-section fixed height wooden frame couch for treatment, Solid Beach Wooden Couch with Adjustable Head Rest, User Weight 225 Kg. Size L 185Cm x W 63 Cm x H 75cm.

ORDERING DETAILS

CODE:	DESCRIPTION
TAB 200 FHW	2 Section Fixed Height Wooden Couch

FIXED HEIGHT WOODEN MAT TABLE



Fixed height Mat Table, Solid Beach Wooden frame, Padding and upholstery fire-resistance rating. Size:200cm x 120cm. 5cm thickness mat, fire retardant. Height:45cm.

ORDERING DETAILS

CODE:	DESCRIPTION
TAB 100 FHSW/M	Fixed Height Mat Table

HEIGHT ADJUSTABLE WOODEN COUCH



Height Adjustable high quality wooden couch with long shelves. Adjustable Head rest and Paper Holder. Length 190cm, Width 76cm, Height: 65 - 84cm.

ORDERING DETAILS

CODE:	DESCRIPTION
TAB 200W/HA	2 Section Height Adjustable Wooden Couch

FIXED HEIGHT COUCH WITH SIDE RAIL



Structure made by Steel Square Tube Frame. Epoxy Powder Coated, 5cm Thick Molded Foam Pad, Fire retardant upholstery, Adjustable Back Rest Assisted by gas spring & Paper Roll Holder with Side Rail. Size : L195cm x W 60cm X H 77cm. Weight Capacity : 200kgs.

ORDERING DETAILS

CODE:	DESCRIPTION
TAB 200 FHS/SR	2 Section Fixed Height Couch with Side Rail