

## TRU TRAC TRACTION UNIT TRACTION COUCH



### ORDERING DETAILS

| CODE:       | DESCRIPTION                                      |
|-------------|--------------------------------------------------|
| <b>4779</b> | <b>TRU TRAC TRACTION UNIT<br/>TRACTION COUCH</b> |

### STANDARD ACCESSORIES

|       |                          |
|-------|--------------------------|
| 48082 | Patient Interrupt Switch |
| 27465 | Patient Data Cards (1x)  |
| 48079 | User Manual on CD        |
| 6870  | TTET-200                 |

### OPTIONAL ACCESSORIES

|       |                     |
|-------|---------------------|
| 66898 | Foot Switch         |
| 86016 | Retractable Casters |

### Economy Fixed Height Traction

|      |          |
|------|----------|
| 6840 | TTFT-200 |
|------|----------|

### SPECIFICATIONS

|                          |                                                                                                                      |
|--------------------------|----------------------------------------------------------------------------------------------------------------------|
| Tru-Trac Traction Device |                                                                                                                      |
| Mains Power:             | 100 - 240 V, 50/60 Hz                                                                                                |
| Power Consumption:       | 75 W                                                                                                                 |
| Current Consumption:     | 3.2 Amps Max                                                                                                         |
| Weight:                  | 14 kg (30 lbs)                                                                                                       |
| Dimensions*:             | 45 x 24 x 45 cm                                                                                                      |
| Electrical Safety Class: | Class 1, Type B                                                                                                      |
| Safety Tests:            | IEC/UL/EN 60601-1, IEC/<br>EN 60601-<br>1-2, CAN C22.2 No.<br>601.1-M90<br>w/A2 - Meets MDD<br>93/42/EEC,<br>CE 0413 |

### FEATURES

- Tru-Trac Traction Unit
- Static, intermittent, and cycling Traction
- User-defined hold, rest and treatment times
- 10 User-defined protocol slots to store customized settings
- Digital monochrome touch screen interface that pivots 270°
- Easily store treatment parameters, pain scales, and pain maps on Patient Data Cards

## IDEAL PLUS TRACTION PACKAGE



### ORDERING DETAILS

| CODE:             | DESCRIPTION                                                                                               |
|-------------------|-----------------------------------------------------------------------------------------------------------|
| <b>TAB 400PKG</b> | <b>IDEAL PLUS TRACTION<br/>(TRU TRAC UNIT WITH HARNESS<br/>FLEXION STOOL AND CERVICAL<br/>ATTACHMENT)</b> |
| <b>TAB 400-B</b>  | <b>IDEAL PLUS TRACTION COUCH</b>                                                                          |

### SPECIFICATIONS

|                     |                |
|---------------------|----------------|
| Minimum height      | 470mm / 18.1"  |
| Maximum height      | 910mm / 35.8"  |
| Cushion width       | 700mm / 27.6"  |
| Head section length | 410mm / 27"    |
| Mid section length  | 680mm / 26.8"  |
| Leg section length  | 740mm / 29.14" |
| Total length        | 2243mm / 88.3" |
| Head tilt angle     | +80" to - 20"  |
| Leg tilt angle      | +70" to 0"     |
| Lift mechanism      | Electric       |
| Breathing hole      | 1              |
| Retractable wheels  | Standard x 4   |
| Weight rating       | 225kg / 496lb  |