

IDEAL PLUS GYNAECOLOGY FIXED COUCH

GP Gynaecology Couch features two electric motors for height and tilt. The backrest is gas spring controlled ensuring easy patient positioning, and enabling the plinth to be moved into Trendelenburg position. Fully adjustable Foot Stirrups are provided as standard, as is the stainless steel operator's tray. This couch can be used as a standard 2 Section Plinth when not being used for Gynaecology.



SPECIFICATION

- Fixed Height Couch with four drawers
- Metal frame and rounded edges for safety
- Adjustable Electric Back and foot tilt
- Dual motor system
- Couch height 80 cm
- Adjustable head section by motor 0 to +80 degree
- Adjustable foot section by Electric Actuator from 0 to +15 degree
- Leg section tilt up to up to 15 degrees.
- Adjustable feet for uneven surfaces
- Soft Material adjustable Knee Crutches
- Removable pillow
- Padding thickness 5cm
- Fire retardant PU Leather

STANDARD FEATURES

- Total Length : 180cm with extension
- Width : 65cm
- Head section length : 88cm
- Middle section length : 65cm
- Leg section length : 50cm
- Height adjustable : Fixed
- Loading weight : 250Kgs
- Upholstery : PVC vinyl with fire resistant
- Certification : CE

ORDERING DETAILS

CODE	DESCRIPTION
TAB 200GFD	Gynaecology Examination Couch with drawer

AVAILABLE COLORS



Manufactured By

