

IDEAL PLUS 300 TRACTION ELECTRIC COUCH /TRACTION PACKAGE



Physiotherapy Traction couch that, features with Rigid Plates are used and reduced the pipes from this new sleek design, twin lift design for ultimate comfort and control. Available in electric. 3 sections Physiotherapy Traction Couch is a fantastic.

STANDARD ACCESSORIES

- Twin lift design.
- Positive/negative head.
- Retractable wheels as standard.
- Breathing hole and plug as standard.
- Electric Hi Lo operation by hand / foot control.
- Sleek design with rigid frame and rounded edges for safety.
- Minimum height 46 cm to maximum height 91 cm.
- Adjustable head section by twin Gas struts.
- Movable foot section during the traction perform.
- Head section has face hole and plug.
- Adjustable feet for uneven surfaces.
- 7.5 cm castors for safety and mobility.
- Retractable dual locking castors on all 4 corners.

TRU TRAC TRACTION UNIT

The Chattanooga Tru-Trac is a simple to use digital touch screen user interface traction unit that offers static, intermittent, and cyclic traction with user definable hold, rest, and treatment times. The Tru -Trac will mount to a platform with minimum dimensions of 300 mm long x 240 mm wide. The clamp can accept a surface/shelf with a maximum thickness of 50 mm.



Code : 7040



Code : 1014



Code : FLEX S400

ORDERING DETAILS

CODE:	DESCRIPTION
TAB300-B	Ideal Plus Traction Couch
TAB 300PKG	Traction Package Inc. Traction Bed, Traction Unit, Traction Attachments, Cervical Traction device with Clevis and Flexion stool.

ORDERING DETAILS

CODE:	DESCRIPTION
4779	Tru Trac Traction Unit