

IDEAL PLUS TILT TABLE



The Ideal Plus variable height tilt table uses dual motors to vary tilt and height adjustment. The very strong and rigid frame fitted with heavy duty locking castors making the tilt table extremely stable at all angle and heights.

FEATURES

- Dual motors operated by hand switch.
- Patient hoist access
- Fitted with tilt angle indicator
- Fitted with quick release mechanism
- Patient table and patient straps as standard

TILT TABLE

The Tilt Table range has been complemented with both 2 Section and Split Leg Models. The 2 Section Tilt Table again has the same functionality as the standard Tilt Table but with an added adjustable backrest, which makes the table extremely versatile, enabling it to be used as a general examination couch as well as a Tilt Table.

FEATURES

- Large Lockable Castors suitable for patient
- Transportation
- Electric Height & Tilt Control
- Patient Handles
- MRSA Resistant Anti Bacterial & Fire
- Retardant Upholstery
- 3 Velcro Harness
- Gas Spring Assisted Back Section
- 180 kg safe working load

ORDERING DETAILS

CODE:	DESCRIPTION
TAB 103E	2 Section Tilt Table
TAB 102E	Split Leg Tilt

SPECIFICATION

• Minimum height	: 790mm / 31.1"
• Maximum height	: 1240mm / 48.8"
• Cushion width	: 650mm / 25.6"
• Total length	: 1960mm / 77.1"
• Head tilt angles	: +87" to 0"
• Lift mechanism	: Electric
• Breathing hole	: No
• Retractable wheels	
• brakes	: Standard x 4
• Weight rating	: 225kg / 496lb
• Net weight	: 88 kg

OPTIONAL EXTRAS

Tab 100 Extra set of patient straps (not illustrated)
"When placing orders please state upholstery colour and optional extras.

ORDERING DETAILS

CODE:	DESCRIPTION
TAB 100	Tilt Table



SPECIFICATION

• Mains Power	: 100 - 240V, 50/60 Hz
• Power	: 75 W
• Current	: 3.2 Amps Max
• Weight	: 14 kg (30 lbs)
• Dimensions*	: 45 x 24 x 45 cm
• Safety Class	: Class 1, Type B

OPTIONAL EXTRAS

Tilting Patient's Table | Tilting Foot Board