

TRU TRAC TRACTION UNIT



FEATURES

- Static, intermittent, and cycling Traction
- User-defined hold, rest and treatment times
- 10 User-defined protocol slots to store customized settings
- Digital monochromatic touch screen interface that pivots 270°
- Easily store treatment parameters, pain scales, and pain maps on Patient Data Cards

SPECIFICATION

▪ Mains Power	: 100 - 240V, 50/60 Hz
▪ Power	: 75 W
▪ Current	: 3.2 Amps Max
▪ Weight	: 14 kg (30 lbs)
▪ Dimensions*	: 45 x 24 x 45 cm
▪ Safety Class	: Class 1, Type B

STANDARD ACCESSORIES

48082	Patient Interrupt Switch
27465	Patient Data Cards (1x)
48079	User Manual on CD

ORDERING DETAILS

CODE:	DESCRIPTION
4779	Tru Trac Traction Unit

TRITON TRACTION UNIT

FEATURES

- Static, intermittent, and cycling traction
- Independent timers allow easy customization of traction pulls in progressive, regressive, and traction phases
- Variable speed motor (30%, 50% and 100%)
- Pretension step eases the patient into traction to prevent muscle guarding and pre-test patient's tolerance to traction force
- Clinical Protocols provided for both cervical and lumbar conditions
- Full color graphic library complete with visual tutorial of belting techniques, anatomical images and common pathologies
- Digital high resolution color touch screen that pivots 270°.
- 80 User defined protocol slots
- Easily store treatment parameters, pain scales, and pain maps on Patient Data Cards and Upgrade with the optional EMG module.



SPECIFICATION

▪ Mains Power	: 100 - 240V, 50/60 Hz
▪ Power	: 75 W
▪ Current	: 3.2 Amps Max
▪ Weight	: 14 kg (30 lbs)
▪ Dimensions*	: 45 x 24 x 45 cm
▪ Safety Class	: Class 1, Type B

STANDARD ACCESSORIES

48083	Patient Interrupt Switch
27465	Patient Data Cards (1x)
48144	User Manual on CD

ORDERING DETAILS

CODE:	DESCRIPTION
4749	Triton Traction Unit