



Auto Traction Table

ZEPU-TQ30

Indications

- Indications for cervical vertebrae: cervical muscle painful spasm, degenerative cervical diseases, cervical disc protrusion, cervical carina nerve root stimulation or compression, degenerative cervical osteoarthritis, vertebrae joint capsule, cervica anterior and posterior longitudinal ligaments
- Indications for lumbar vertebrae: lumbar disc herniation, degenerative lumbar disc diseases, lumbar small joint dysfunction, lumbar degenerative osteoarthritis, lumbar muscle painful spasm and tension.

Features

- Three tracking modes meeting the requirements of treatment
- Primary and secondary tracking forces are flexible
- Offering double protections: internal protection over stretching program and external emergency stop remote control
- Cervical and lumbar vertebrae could be operated at the same time
- Heated bed mattress is one of options

Specification

Weight:	main frame: 15kg; bed: 122kg
Size (Main frame) :	410 x 200 x 480mm
Size (Bed) :	1800 x 1000 x 680-1050mm
Treatment Modes:	5 internal prescriptions and 99 customizing plans
Max Tracking Force:	990Kg±20%
Treatment Time:	0.99±30s
Ascending Speed:	0.9mm/sec-1.8mm/sec
Tolerance:	±10%
Descent speed:	0.9mm/sec-1.8mm/sec

ZEPU-TQ30 /SET

ZEPU-TQ30/SET

ZEPU-TQ30 Electric Traction Bed Set

THE MULTI-POSITION BED COMBINED WITH INTELLIGENT CONTROL SYSTEM MAKES TRACTION THERAPY MORE INTELLIGENT, COMFORTABLE, AND EFFICIENT.

ZEPU-TQ30 provides an intelligent control system, it adopts LCD touch screen with 5 build-in treatment and 99 customizable prescriptions. Multi-position bed configuration, the bed can be smoothly lifted and lowered for traction angle setting.

During the traction process, the bed separates itself to reduce friction loss, and the bed is adjusted in 6-fold sections to provide the best traction position for patients.

Innovative V-shaped occipital design neck traction device makes neck traction safer and more comfortable, and traction position more precise.

