


IDEAL PLUS 2 SECTION WIDE NEUROLOGY PLINTH



The classic dual life wide neurology couch is probably the most stable and rigid couch currently available. Featuring a large upholstered top and an exceptional low operation height, 46cm (18") yet still allowing clearance for mobile hoist procedures.

STANDARD FEATURES

- Positive/negative head section
- Retractable wheels as standard
- Breathing hole and plug as standard
- Electric Hi Lo operation by hand / foot control
- Twin lift design with sleek frame and rounded edges for safety
- Minimum height 46cm to maximum height 91cm
- Adjustable head section by dual gas struts from -20° to +80°
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm castors for safety and mobility
- Padding thickness 5cm



SPECIFICATION

• Total Length	: 205cm
• Width	: 100cm
• Head section length	: 69cm
• Leg section length	: 136cm
• Height adjustable	: 46cm to 91cm
• Loading weight	: 225Kgs
• Upholstery	: PVC vinyl with fire resistant
• Certification	: CE

FUNCTIONS

- Adjustable backrest (+ 65°)
- Includes the following accessories:
- Breather hole
- 4 large, retractable wheels

ORDERING DETAILS

CODE:	DESCRIPTION
TAB 200W / ELEC	2 Section Wide Neurology Plinth

AVAILABLE COLORS

Black	Grey	Beige	Vivid Blue	Spruce Blue	Navy Blue