

PORTABLE IONTOPHORESIS

The Pocket Physio Line includes a new electrostimulator specific for analgesic and anti-inflammatory applications: the IONOTENS model. Its software includes 33 pre-set programs for the most common applications and 4 free programs for setting personalised programs. It issues analgesic tens, ionophoresis, Iontoforesis and anti-inflammatory currents.

SPECIFICATION

- Power Supply : Rechargeable battery / power mains
- No. Output Channel : 1
- No. of Program (Hz) : 33 pre-set + 4 free Programs
- Emission : Continuous, modulated
- Frequency : 1-20 (HZ)
- Dimensions : 18 x 11 x 4cm, 400gr only device,



ORDERING DETAILS

CODE:	DESCRIPTION
BIO-1026	Portable Iontophoresis Therapy

PORTABLE MAGNETOTHERAPY

New Pocket Magneter is the new low-frequency magnetotherapy device, easy to use with reduced dimensions. New Pocket Magneter is ideal for medical and home applications with its 68 pre-set programs for the treatment of the most common pathologies. New Pocket Magneter also contains 20 free programs to realize personalized therapies for all kind of patients.

SPECIFICATION

- Power Supply : Rechargeable battery / power mains
- No. Output Channel : 2
- No. of Program (Hz) : 68 pre-set + 20 free Programs
- Emission : Continuous, modulated
- Frequency : 120 gauss (60 Gauss for each solenoid).



ORDERING DETAILS

CODE:	DESCRIPTION
BIO-1011	Portable Magnetotherapy

PORTABLE PHYSIO URO

The Pocket Physio Line includes a new electrostimulator specific for urogynecological applications: the URO model. For the pelvic functional recovery the device is endowed with special probe for personal use. Its software includes 20 pre-set programs for the most common applications and 4 free programs for setting personalised programs. It is so easy to use that the patient can make applications alone at home (2 waveforms) (Vaginal Probe).

SPECIFICATION

- Power Supply : Rechargeable battery / power mains
- No. Output Channel : 1
- No. of Program (Hz) : 33 pre-set + 4 free Programs



- Frequency : 1-20 (HZ)
- Dimensions : 18 x 11 x 4cm, 400gr only device,
- Weight : 1.5 kg with accessories

ORDERING DETAILS

CODE:	DESCRIPTION
BIO-1058	Portable Physio Uro With V. Probe