

MOBILE 2 RPW SHOCKWAVE

FEATURES

The Mobile 2 RPW is the prodigal son of the Mobile RPW, with improved versatility (a broader range of shock frequency) and greater power (up to 5 bar of pressure). These slight, yet significant, adjustments make it an all round better tool for helping patients recover, treating resilient conditions and more effectively administering stimulation to injured tissues. The versatility and usefulness of the Mobile 2 RPW is such that it is recommended for all of the following environments and applications: orthopaedics, physiotherapy, aesthetic therapy, chiropractors, rheumatology, sports clubs, veterinary, physical medicine, acupuncture and podiatrists.

- All the features of the original Mobile RPW but with improved functionality.
- Single and Continuous shocks: Frequency
- D-actor: 1 - 21 Hz Frequency V-actor: 1 - 31 Hz.
- 1-5 Bars of pressure, infinitely variable with precise adjustments of 0.1 increments.

STANDARD ACCESSORIES

- 32000 R-SW Handpiece Set, including: R-SW Hand Piece, Projectile (x2), Guiding Tube (x2) Allen Key, Sealing Set, Wrench, Cleaning Brush, Case for R-SW Hand Piece
- 23825 Revision Kit, Including: Projectile (x2), Guiding Tube (x2), Sealing Set (x2), O-Ring Guide, Case for Revision Kit
- 28178 Ro40 'Beam' transmitter
- 29724 D20-S R-SW 20 mm transmitter
- 23825 Revision kit
- 13-5182 Conductor transmission gel 250 ml

OPTIONAL ACCESSORIES

- 22376 Shaft cushion
- 22375 Handle cushion
- 28178 Ro40 'Beam' transmitter
- 29801 A6 shock transmitter
- 18638 R15 shock transmitter
- 29729 C15 shock transmitter
- 29726 F15 shock transmitter
- 29728 DI15 shock transmitter
- 29724 D20-S shock transmitter
- 28736 D20-T shock transmitter
- 29725 D35-S shock transmitter



- 29802 T10 shock transmitter
- 28730 V-ACTOR II Handpiece Set (See p 23)
- 28740 V25 V-ACTOR® 25 mm Vibration Transmitter
- 28741 V40 V-ACTOR® 40 mm Vibration Transmitter
- 2910 Mobile 2 RPW Shockwave

SPECIFICATION

- R-SW Operating Mode : Single shock, continuous shock 1.0 - 4.0 bar / max. 21 Hz 4.6 bar / max. 15 Hz 5.0 bar / max. 12 Hz
- V-ACTOR operating mode : 1.4 - 3.0 bar / max. 31 Hz 3.8 bar / max. 22 Hz 5.0 bar / max. 13 Hz
- Mains Power : 100 - 240 VAC
- Mains Frequency : 50 - 60 Hz
- Mains Fuse : T4AL/250 VAC
- Power Consumption : Max, 200 VA
- Compressed Air Output : 1-5 bar
- Ambient Temp. During Operation : 10-40 °C
- Ambient Temp. Storage & Transport: 0- 60 °C Frost-free
- Ambient Air Pressure: 800-1,060 h Pa
- Classification According to MDD: Class 2A Device
- Ingress Protection Rating: IPX1
- Weight Base unit: 8.2 kg
- Dimensions Base unit: 32 x 32 x 21 cm

ORDERING DETAILS

| CODE:    | DESCRIPTION                     |
|----------|---------------------------------|
| CHA 2910 | Mobile 2 Portable RPW Shockwave |