

INTELECT® MOBILE LASER



FEATURES

- Class 3B Low Level Laser Therapy
- Independent control over all parameters
- Continuous and pulsed operation available
- Fully variable pulsed frequencies
- Real time feedback of energy delivered
- Display of energy in choice of joules or joules/cm²
- Helpful menu of Clinical Indications specific to the laser applicator
- 15 user-defined memory slots on board
- Built in acupuncture point finder
- 17 different applicators available, including LED, laser diodes, and clusters
- Innovative design allows for cart, table top, wall or mobile use
- Backlit LCD display for better visibility in low light settings
- Improved pq button for easier navigation

LASER APPLICATORS (Singles)

27799	670nm LED Diode 10mW
27803	850nm Laser Diode 40mW
27840	850nm Laser Diode 100mW
27804	850nm Laser Diode 150mW
27841	850nm Laser Diode 200mW
27805	820nm Laser Diode 300mW

9 Diode Cluster

27810	9 Diode Cluster Applicator, 290mW Total: 5 x 850nm 50mW Lasers, 4 x 670nm 10mW LED
27812	9 Diode Cluster Applicator, 1040mW Total: 5 x 850nm 200mW Lasers, 4 x 670nm 10mW LED

13 Diode Cluster

27813	13 Diode Cluster Applicator, 265mW Total: 3 x 850nm 50mW Lasers, 7 x 670nm 10mW LED, 3 x 950nm 15mW LED
27814	13 Diode Cluster Applicator, 415mW Total: 3 x 850nm 100mW Lasers, 7 x 670nm 10mW LED, 3 x 950nm 15mW LED
27816	13 Diode Cluster Applicator, 715mW Total: 3 x 850nm 200mW Lasers, 7 x 670nm 10mW LED, 3 x 950nm 15mW LED

19 Diode Cluster

27815	19 Diode Cluster Applicator, 325mW Total: LED only, 6 x 670nm 10mW LED, 7 x 880nm 25mW LED, 6 x 950nm 15mW LED
-------	--

33 Diode Cluster

27809	33 Diode Cluster Applicator, 565mW Total: LED only, 12 x 670nm 10mW LED, 13 x 880nm 25mW LED, 8 x 950nm 15mW LED
27802	33 Diode Cluster Applicator, 690mW Total: 5 x 850nm 50mW Lasers, 12 x 670nm 10mW LED, 8 x 880nm 25mW LED, 8 x 950nm 15mW LED
27807	33 Diode Cluster Applicator, 940mW Total: 5 x 850nm 100mW Lasers, 12 x 670nm 10mW LED, 8 x 880nm 25mW LED, 8 x 950nm 15mW LED
27808	33 Diode Cluster Applicator, 1440mW Total: 5 x 850nm 200mW Lasers, 12 x 670nm 10mW LED, 8 x 880nm 25mW LED, 8 x 950nm 15mW LED



ORDERING DETAILS

CODE:	DESCRIPTION
2779	INTELECT® MOBILE LASER

STANDARD ACCESSORIES

27525	Laser Protection Glasses (2x)
28063	User Manual on CD

OPTIONAL ACCESSORIES

27904K	Remote Interlock Kit
27467	Intelect Mobile Carrying Bag
27478	Battery Pack
2775	Therapy System Cart
2883	Cart Adapter

SPECIFICATIONS

Main Power:	120-240 VAC, 50 60 Hz
Weight:	2.3 kg
Dimensions*:	33 x 29 x 16.3 cm
Electrical Safety Class:	Class 1, Type BF
Safety Tests:	IEC 60601-1, IEC 60825-1, EN 60601- 2-22

