

Cold Therapy





Cryo/Cuff™ IC Cooler

Indications

- Provides intermittent compression to minimize hemarthrosis, swelling and pain

Features and Benefits

- Integrated pneumatic pump within the cooler lid provides a 45 second on/30 second off cycle of automated compression and cold therapy
- Its ease of operation makes it ideal for post-operative recovery, trauma, athletic training rooms and home use
- The IC Cooler is compatible with all Cryo/Cuff wraps

Ordering Information

Cryo/Cuff IC Cooler	
Part Number	Description
52A	Cryo/Cuff IC Cooler

Accessories	
Part Number	Description
25-0238	Cryo/Cuff IC Lid
25-4882	Cryo/Cuff IC Cooler power supply (US)
25-2740	Universal power supply (all markets excl. US)



Cryo/Cuff™ Cooler

Indications

- General haemarthrosis
- Swelling and pain post surgery and post injury

Features and Benefits

- Cryo/Cuff utilizes gravity to deliver chilled water to anatomically specific cuff
- Cooler holds water and ice needed for 6 to 8 hours of cryotherapy
- Compatible with all Cryo/Cuff pads
- Includes tube assembly and insulation disk
- Labelled with simple-to follow directions for use

Ordering Information

Cryo/Cuff Cooler	
Part Number	Description
2125	Cryo/Cuff Cooler

AutoChill® System

Indications

- Designed for use with the Aircast Cryo/Cuff Cooler
- During operation, the electronic pump automatically cycles and compresses the air in the cooler to force chilled water into the cuff

Features and Benefits

- Automatically cycles the water between the cooler and cuff
- Repetitive cycling maintains continuous cold and applies soothing pulsating pressure
- Eliminates the need for manual recycling

Ordering Information

AutoChill System	
Part Number	Description
20B	AutoChill System (120V AC) - US Power Supply
20BG-E	AutoChill System (230V AC) - European Power Supply
2007	Replacement tube 50 cm
2007XL	Replacement tube 1.50 m



Cryo/Cuff™ Wraps & Cuffs

Indications

- The anatomically designed cuffs provide the therapeutic benefits of controlled cold compression to minimise haemarthrosis and swelling and reduce pain

Features and Benefits

- Cuffs are anatomically designed to cover all major joints of the body
- Cuffs provide maximum coverage for optimum cryotherapy treatment
- Soothing cold covers the affected area reducing swelling and providing pain relief
- Cuffs are simple to apply and remove and may be detached from the cooler at any time for mobility without interruption of treatment



Ordering Information

Cryo/Cuff Wraps & Cuffs					
Part Number		Description	Size		
Cuff	Cuff and Cooler				
10C01	10C	Foot Cryo/Cuff	M	Foot Circumference	23-33 cm (9"-13")
10B01	10B	Foot Cryo/Cuff	L		25-43 cm (10"-17")
10A01	10A	Ankle Cryo/Cuff	U	Universal	
11C01	11C	Knee Cryo/Cuff	S	Circumference 15 cm (6") above patella	25-50 cm (10"-19")
11A01	11A	Knee Cryo/Cuff	M		45-60 cm (18"-23")
11B01	11B	Knee Cryo/Cuff	L		50-80 cm (20"-31")
11ASC	-	Knee Cryo/Cuff SC	U	Universal	
13C01	-	Calf Cryo/Cuff	U	Calf Circumference	36-51 cm (14"-20")
13A01	-	Thigh Cryo/Cuff	L	Thigh Circumference	43-58 cm (17"-23")
13B01	-	Thigh Cryo/Cuff	XL		51-69 cm (20"-27")

Cryo/Cuff Wraps & Cuffs					
Part Number		Description	Size		
Cuff	Cuff and Cooler				
16A01	16A	Hand/Wrist Cryo/Cuff	Universal		
15A01	15A	Elbow Cryo/Cuff	Universal		
12A01	12A	Shoulder Cryo/Cuff	U	Torso Circumference	81-122 cm (32"-48")
12B01	-	Small Shoulder Cryo/Cuff	S		< 86 cm (< 34")
12AXL01	12AXL	Shoulder Cryo/Cuff with Extra Long Strap	XL		107-137 cm (42"-54")
1220XL	-	Shoulder Strap	XL		
14A01	14A	Back/Hip/Rib Cryo/Cuff	Universal		
11P01	-	Pediatric Knee/ Elbow Cryo/Cuff	Pediatric		
10P01	10P	Pediatric Ankle Cryo/Cuff	Pediatric		