

QMD® CRYO THERMAL PT

qmd enlarges the appliance family! After the success[®] of the studio version, thanks to its innovations and ease of use, we decided to satisfy even more the needs of our customers. We introduced the portable version of our cryo-thermal: qmd cryo thermal PT.

The qmd cryo thermal PT is a portable therapy system for the application of cryotherapy, thermotherapy, thermal shock and the exclusive contrast therapy. With an ergonomic handpiece that allows direct skin contact, the desired therapeutic effects can be specifically obtained.

DESIGN

qmd cryo thermal PT does not betray the qmd philosophy and always maintains its unmistakable design, which is characterized by simple, essential

FUNCTIONALITY

We made a device with really small dimensions. It can be easily transported, even outside your studio and is able to meet any expectation.

MODULARITY

We wanted to maintain modularity as one of the important features of our brand: qmd[®] cryo thermal PT can be upgraded at a later stage with the qmd[®] ultrasound module.

Even in the portable version, the applicable therapies are not reduced:

- cryo therapy (-2 °C up to +15 °C) ❄️
- thermo therapy (+30 °C up to +42 °C) 🔥
- contrast therapy
(alternation of heat and cold) ❄️ 🔥 ❄️ 🔥
- thermal shock ❄️ 🔥

TECHNICAL DATAS

- Dimensions 24 x 34 x 22 h cm
- Weight 7.5 kg
- 10.1" display
- Water cooling system
- USB input to update your software
- Ethernet input for remote control
- No condensation water, no gas and no jet of air

Pathologies treatable with the cryo thermal PT are:

- Rheumatoid arthritis and inflammatory
- Diseases in general
- Back pain
- Osteoarthritis
- Muscle problems
- Acute injuries
- Post-operative convalescence
- Exploiting the effects of heat and cold it is possible to treat many other pathologies.

**ORDERING DETAILS**

CODE:	DESCRIPTION
Q-0203	Cryo Thermal Therapy