

Digital Air Compression Therapy System

Perfect for pre and post workouts, as well as rehabilitation, the sequential compression pump effectively increases circulation to reduce swelling and leg pain, venous edema, lymphedema, post-thrombotic syndrome, peripheral arterial disease, diabetes and other conditions. Gentle air compression therapy improves blood circulation and enhances the pressure of muscular tissues to relieve sore muscles.



Operation panel advantages:

- The ability to turn off individual chambers in the cuffs, e.g. to avoid injured places
- Possibility to change the parameters of the treatment while the device is operating
- Automatic calculation of gradient values after setting minimum and maximum values
- Available for the visually impaired, set voice prompts, customizable special keyboard
- System can be updated and upgraded for continuous access to the latest features
- Allows immediate termination of the treatment and deflation of air from the cuffs
- Real pressure measurement during the treatment
- Adopt touchable screen, easy to operate and provide a clear display
- Digital Air Compression Therapy System
- Treating time: 10-30mins/ always on
- Power: 30W/65VA
- Pressure range : 30-250 mm Hg
- Mode : 6 modes
- Dvice size : 24cm(L)*20.5cm(W)*20cm(H)
- Unit including 1pair full leg sleeves. Waist and Arm Sleeves are Option



CODE: F 3001D

DESCRIPTION: DIGITAL AIR COMPRESSION THERAPY SYSTEM