

DIGITAL TROPHIC STIM



ORDERING DETAILS

CODE	DESCRIPTION
Em-2100	Trophic Stimulator W/24 Programmes

The STIM 8 Trophic Stimulator is an advanced micro-processor controlled two channel combination unit of TENS(Transcutaneous Electrical Nerve Stimulator) and EMS(Electrical Muscle Stimulator) offering 8 trophic stimulator treatments.

Program	Intended for Use
1-3	Facial Palsy
4-6	Facial Synkinesis
7-9	Brachial Plexus Injury
10-12	Trophic Stimulation for Stroke Patients
13-15	Arthritic Conditions
16-18	Knee
9-21	Pain Management
22-24	Sports & Recent Injury

DIGITAL OBSTETRICS TENS



ORDERING DETAILS

CODE	DESCRIPTION
EV-508P	Portable Digital Obstetric Tens

Birth Assist Obstetric TENS provides with another option, encouraging your body to produce a greater amount of its natural painkillers (endorphins). TENS has no known side effects and is a proven method of managing pain during labour. This obstetric TENS machine is a dual channel unit that is both durable and easy to use. .

- Dual channel unit
- Output can be adjusted to each pair of electrodes, allowing for response to degree of local pain in the user
- Burst mode – allows extra pain relief
- Burst mode is ideal for pain of contractions
- Includes hand held trigger for burst mode
- Digital display for clear, accurate and easy use
- User friendly design for simple operation
- Natural alternative to pain relieving drugs

DIGITAL INCONTINENCE STIM



ORDERING DETAILS

CODE	DESCRIPTION
2400A	Portable Digital Incontinence

The EM-2400 incontinence stimulator is a battery operated pulse generator that stimulates the muscles of the pelvic floor by means of a vaginal or rectal electrode.

Channel	Single channel
Pulse Amplitude	Adjustable, 0-80mA, Max output 80mA (peak to peak) into 500ohm load.
Output Voltage	Adjustable, 0-40V, Max output 40V (peak to peak) into 500ohm load.
Wave Form	Asymmetrical Rectangular Biphasic Pulse.
Power Supply	Two AA Batteries
Pulse Rate	Pre-setting, 1-100 Hz
Pulse Width	Pre-setting, 150 μ s, 200 μ s or 250 μ s
Program	A, B, C, D, E, five pre-set program of options.