


EXPIA EXPERT - 2 SECTION EXAMINATION COUCH



Expia Expert - treatment tables which combine simplicity, quality and design. The shape of the base is inspired on a runner who is down on his marks, just about to start the race; ready and steady to go. Tables conceived on the premise that the best ideas are the simplest ones. Electric 2-section therapy table for treatment on patients lying down in supine or prone positions or semi-sitting position. The backrest is adjustable from 0° to 65°. The table is equipped with 1 motor that controls the table height through the perimeter foot switch placed at the base of the table. Electric examination bed with electric height adjustment from 47 - 85 cm. Operated by remote control or foot switch. Removable Face Cushion.



SPECIFICATION

▪ Motor	: 1
▪ Dimensions	: 196 x 68 cm including hand rest
▪ Height	: 47-85 cm
▪ Weight	: 95 kg
▪ Thickness	: 5cm
▪ Density	: 40
▪ Load Capacity	: 200 kg
▪ Power supply	: 230V - 50/60 Hz

FUNCTIONS

- Adjustable backrest (0/+ 85°)
- Lifting mechanism from both front and rear
- 2 wheels for easy movement
- Security lock

ORDERING DETAILS

CODE:	DESCRIPTION
EXP 220	Expia Expert 2 Section Examination Couch