

EXPIA EXPERT - 3 SECTION TREATMENT COUCH



Expia Expert is a collection of treatment tables which combine simplicity, quality and design. The shape of the base is inspired on a runner who is down on his marks, just about to start the race; ready and steady to go. Tables conceived on the premise that the best ideas are the simplest ones. Their soft, round and fine shapes contrast their strong personality and born-fighter spirit.



SPECIFICATION

• Product size	: 200 x 68 cm
• Height	: 54 / 95 cm
• Motor	: 1nos
• Weight	: 125 kg
• Thickness	: 4 cm
• Density	: 40
• Breather hole	: yes
• Paper roll holder	: yes
• 2 wheels turning	: 360°
• Loading capacity	: 250 kg
• Operation	: Peripheral swiss Linak High Reliability' actuators
• Size	: 1900 x 6350 x 540 - 960mm

OPTIONAL

Heating, Safety lock , Bed cover , 2 Motors Height + Seat

ORDERING DETAILS

CODE:	DESCRIPTION
EXP 303/1M	Expia 3 Section Couch

ORION EXAMINATION COUCH WITH 3 MOTORS



3-section physiotherapy table intended for massage, functional rehabilitation and osteopathy treatments.

OPTIONAL

Heating, Safety lock , Bed cover , Handlebar Support

SPECIFICATION

• Product size	: 198 x 68 cm
• Height	: 54 / 95 cm
• Motor	: Safety Clutch, Safety button with key (optional)
• Weight	: 125 - 135 kg
• Thickness	: 5 cm
• Density	: 40
• Breather hole	: Yes
• Paper roll holder	: Yes
• 4 wheels turning	: 360°
• Loading capacity	: 250 kg
• Operation	: Peripheral swiss Linak High Reliability' actuators
• Size	: 1900 x 680 x 540 - 950mm

ORDERING DETAILS

CODE:	DESCRIPTION
ORN 301/2M	Expia 3 Section Couch