

MTS Z100

Magnetic Stimulation System

-  Incontinence
-  Sexual Dysfunction
-  Chronic Pain
-  Pelvic Organ Prolapse



1st-line Non-invasive Solution

Pictures are for reference only. The actual product may be different from region to region.


OVERVIEW

MTS Z100 is a **pelvic magnetic stimulation** (PMS) system that applies **targeted and powerful** pulsed magnetic field on deep pelvic muscles and sacral nerves to treat urinary & fecal incontinence, sexual dysfunction, chronic pain, and mild pelvic organ prolapse for **men, women, and children**. It is a completely **non-invasive** option that allows patients to stay **fully clothed** throughout the treatment.

MEN & WOMEN & CHILDREN

FULLY CLOTHED


FREE OF OPERATION



TECHNOLOGY

- ◀ **FMS:** Functional Magnetic Stimulation, the system generates pulsed magnetic field to cause pelvic floor muscles to performance thousands of effective contractions without direct contact with electrodes or probes.
- ◀ **SNMS:** Sacral Nerve Magnetic Stimulation, the system sends mild electromagnetic stimulations on sacral roots to regulate pain and stress, improving urinary and bowel control.

ADVANTAGES



Outstanding Performance

**7 Tesla
magnetic field**
for
deep penetration

**100Hz
pulse frequency**
ensures
effective treatment

**Liquid
cooling system**
for non-
intermittent work

All Genders

- Benefits for Male, Female, and Children

Adjustable Positions

- Allow effective stimulation in different areas

MEP Test

- Guide proper stimulation target

Convenient & Efficient

- Patients remain fully dressed, free of operation
- Predefined & customizable session plans for each patient
- Comfortable arm-chair design enables pleasant treatment experience



Pictures are for reference only. The actual product may be different from region to region.