

TESLA MAGNETO THERAPY



Tesla is an ELF magneto therapeutic unit designed for professional applications. Tesla generates magnetic field(ELF) characterized by: Low frequency, Negligible electric field, Low radiated power, Negligible thermal effects, High tissue penetration. The unit has a graphic touchscreen display which allows the interactive control of the unit. Tesla 2 can have at same time connected a couple of plain transducers and a cylindrical transducer or two cylindrical transducers. The outputs can be activated either contemporarily or one a time.

FEATURES

- Output frequency 2 – 100 Hz
- Max output level 100 / 10 Gauss/mT
- Power supply 90 - 240 Vac / 50-60 Hz
- Electrical input power
- 120 VA / 150 VA

ORDERING DETAILS

CODE:	DESCRIPTION
TSL001MG	Tesla Magnetotherapy
TSL001WB	Magneto Unit with Wooden Bed & 60cm Cylinder

TESLA PLUS MAGNETO THERAPY

Tesla pulse is a device designed for professional applications for the treatment of oste-articular pathologies. Tesla pulse can command a pair of flat transducers (optional) and a new handle of pulsed magnetic field. Outputs can be activated individually. Tesla Pulse in PULSE mode through the interchangeable circular head of the handpiece emits a pulsed magnetic field characterized by: Low frequency, negligible electric field, high magnetic pulsed field (1 Tesla) (Radiated Magnetic Pulse RAMP), negligible thermale effects, high penetration in tissues, this particular field results in non-invasive and non-painful applications, favoring draining and cellular normalization through tissue and cell biostimulation for the treatment important traumatic pathologies with or without edema, inflammatory and or painful forms such as cervical arthrosis, lumbagoe, scapolomeric periarthritis, sciatica, tendinitis, sports trauma (stretching, tearing, etc.) Tesla Pulse has a touchscreen graphic display that allows interactive control of the unit.

FEATURES

- Power supply 90 - 240 Vac / 50-60 Hz
- Electrical input power 120 VA
- FEATURES pulse mode
- Output frequency 0.5– 5Hz
- Max output level 10 kGauss (1T)



ORDERING DETAILS

CODE:	DESCRIPTION
LTSL001PMG	Tesla Magnetotherapy