

INTELECT MOBILE 2 STIM



Intellect® Mobile 2 Stim is a 2-channel electrotherapy system. It has 5 types of current: interference, VMS, high voltage, asymmetrical, biphasic and symmetrical biphasic TENS current. The different current types can be quickly selected via the simple user guidance and touch screen. Instructions for the respective electrode placement are also displayed for the selected treatment area.

STANDARD FEATURES

- Full color 7" high resolution touch screen
- Upgradable software : future proof technology
- Import / Export capability (Bluetooth)
- Treatment data storage
- Lightweight and ergonomically designed handle
- Home use: compliant with IEC 60601-1-11

KEY FEATURES - INTERFACE

Intelligent interface with touch buttons and interactive lights indicating status of device or port. Suggested Protocol Setting (SPS): built in suggested protocols based on currently available clinical evidence. Anatomical library with high resolution images for patient explanation. Custom protocols can be stored

WAVEFORMS

- Interferential current (Traditional IFC, 4-pole)
- Premodulated IFC (2-pole)
- Asymmetrical & Symmetrical Biphasic current - TENS
- High Voltage Pulsed Current (HVPC)
- Isoplanar, Dipole
- **VMS, VMS Burst & VMS FR**
- Monophasic Triangular & Rectangular Pulsed
- Galvanic & Galvanic Interrupted
- Low Level Galvanic
- Russian & Träbert Current
- Han Stimulation & Microcurrent
- Diadynamic LP, MF, DF, Cpid/Cpiso
- Surged Monophasic Triangular & Rectangular



STANDARD ACCESSORIES

- Rubber electrode 6 x 8 cm (1/pack)
- Stim Ch 1/2 Leadwire Kit
- Dura-Stick® Plus 5 cm (2") cloth, square
- USB Drive Intellect Mobile 2
- Printed Quick Start Guide Mobile 2
- Chattanooga Straps
- Power cord EU Intellect Mobile 2

OPTIONAL ACCESSORIES

- Mobile 2 cart
- Battery Pack
- Stim 1/2 lead wire kit
- Stim 1/2 XL lead wire kit
- Hi volt Probe
- Carry Case
- Vaccum Module

SPECIFICATION

• Main input voltages	: 100-240V AC, 1.0 to 0.42A, 50/60Hz
• Main input voltages	: 100-240V AC, 1.0 to 0.42A, 50/60Hz
• Electrical Class	: Class II
• Electrotherapy	: Type BF
• Mode of operation	: Continuous
• Dimensions / Weight	: 340 x 355 x 150 mm / 3.1Kg
• Conformity	: This device complies with the Applicable standards EN/IEC 60601-1 & 60601-1-11
• CE 2797	: Classification according to MDD class IIa device 93/42 60601-1

ORDERING DETAILS

CODE:	DESCRIPTION
CHA 15-1203	Intellect Mobile 2 Stim
CHA 15-1210	Mobile 2 Vaccum Module with