


ORION - OSTEOPATHY 2 SECTION EXAMINATION COUCH



2 section therapy table with hand rest. 1 motor that controls the height of the table through the built-in perimeter foot switch. Equipped with Twin Lift: the patient can adjust the armrests to reduce tension and get shoulder to relaxed during treatments. 4 wheels turning 360°.

Peripheral Switch

Ergonomic and silent height-controlling device all around the table. The mechanism in each side consist of two bars fully integrated in the base. In the front and rear part of the table, the height motion is operated by means of two foot switches for both up and down movements respectively.



SPECIFICATION

▪ Motor	: 1 Engine (height)
▪ Dimensions (lxb)	: 200 x 55-80 cm including hand rest
▪ Height (Min. & Max.)	: 55/95 cm
▪ Weight	: 141 kg
▪ Thickness	: 5cm
▪ Density	: 40
▪ Load Capacity	: 250 kg

TECHNICAL DATA

- Breather hole
- Handlebar support
- Wheels 360°
- Peripheral switch
- Security lock

ORDERING DETAILS

CODE:	DESCRIPTION
ORN 201/OSP	Orion Osteopathy Couch