

## Rehab



| Ordering Information                             |                                                                                 |
|--------------------------------------------------|---------------------------------------------------------------------------------|
| Part Number                                      | Description                                                                     |
| Rehab (Full kit includes 4 bags of 4 electrodes) |                                                                                 |
| 2533110                                          | Rehab EN-FR-IT-DE-SP-NL Plug EU                                                 |
| 2533111                                          | Rehab EN-TU-PT-NL-SP Plug EU                                                    |
| 2533112                                          | Rehab EN-DA-SV-NO-FI Plug EU                                                    |
| 2533113                                          | Rehab EN-CZE-PL-RU Plug EU                                                      |
| 2533114                                          | Rehab EN-FR-IT-DE-SP-NL Plug US                                                 |
| 2533115                                          | Rehab EN-FR-IT-DE-SP-NL Plug UK                                                 |
| 2533117                                          | Rehab EN-FR-IT-DE-SP-NL Plug AUS                                                |
| Standard Accessories                             |                                                                                 |
| 683010                                           | Fast charger                                                                    |
| 941213                                           | Battery pack                                                                    |
| 601132                                           | Set of 4 pin cables                                                             |
| 6111944                                          | Snap adaptor kit                                                                |
| 690001                                           | Protecting cover                                                                |
| 45264XX                                          | User manual and practical guide on CD/USB                                       |
| 885932                                           | Quick start guide + warning leaflet                                             |
| 4526645                                          | Electrode placement leaflet                                                     |
| 602047                                           | Bottle of gel                                                                   |
| 980020                                           | Motor point pen                                                                 |
| 680055                                           | Transportation pouch                                                            |
| 949000                                           | Belt clip                                                                       |
| Optional Accessories                             |                                                                                 |
| 6601160                                          | Mi-sensor lead wire (1pc)                                                       |
| 690001                                           | Antishock protection sleeve                                                     |
| 111.970                                          | External battery charger                                                        |
| Technical Specifications                         |                                                                                 |
| Channels:                                        | 4                                                                               |
| Constant current:                                | Up to a resistance of Max 1300 Ohm/ channel                                     |
| Maximum pulse intensity:                         | 120 mA                                                                          |
| Pulse intensity increments:                      | 0.25 mA                                                                         |
| Pulse width:                                     | 30 to 400 $\mu$ s                                                               |
| Pulse frequency:                                 | 1 to 150 Hz                                                                     |
| MDD:                                             | Class IIa                                                                       |
| IP class:                                        | IP22                                                                            |
| Applied part:                                    | Type BF<br>Comply with safety tests: IEC 60601-1, IEC 60601-1-2, IEC 60601-2-10 |
| Battery 941213:                                  | Nickel metal hydride (NiMH) 4.8 V 2000 mA/h                                     |
| Dimensions:                                      | 138 x 95 x 33mm                                                                 |
| Weight:                                          | 300g battery included                                                           |



### A Standard for Rehabilitation

REHAB is a 4 channel stimulator developed for Physiotherapists and Rehab Clinics that require a standard wireless option of electro stimulation. It meets all of the requirements a professional demands; pain management, treatment of vascular diseases, and muscle stimulation with basic programs.

### Features

#### Programs

Four standard program categories: Rehabilitation I, Pain Relief I, Vascular and Conditioning I. Standard program configuration.

#### Muscle Intelligence™ Technology

- Mi-Ready (integrated Mi Technology on all 4 channels, but Mi-Cable not included)
- 4 functions: Mi- Scan, Mi-Tens, Mi-Range, Mi-Action

#### 2+2

Choose one program for channel 1 and 2, and then select another program from the 2+2 program list for channel 3 and 4. This allows treatment of two different body parts, or two patients simultaneously, or combines NMES and TENS as one treatment.

#### Skip

Manage the session time and move directly to the next phase (recovery) when needed.

#### Motor Point Pen

This very useful tool in combination with the dedicated program helps to ensure the optimal electrode placement.

#### Programmability

Creates and personalizes programs in order to customize the patient specific parameters.

#### Detailed Practical Guide

Guidelines and treatment recommendations as well as technical parameters for all programs included in the stimulator are provided. The fold out cover includes illustrations of the suggested electrode placements.

#### Autonomy and Speed Recharge

Daily use autonomy with the charger included in the kit, a full cycle of recharging will take a maximum of 2.5 hours.

#### Warranty

The warranty is for 1 years from the date of purchase.

#### Waveforms



- Symmetrical Biphasic 100% compensated