


SONORA COMBINATION



SONORA COMBINATION 2 CHANNEL ELECTROTHERAPY AND 1 & 3 MHz ULTRASOUND WITH 7" TOUCH SCREEN.

SONORA COMBI with automatic tuning to use ultrasound transducers in the range 1 MHz – 3 MHz together with

- Electrotherapeutic current.
- Through Sonora it is possible to obtain delivery of continuous or modulated ultrasound with programmable intensity and delivery time to fit the therapeutic deep treatment for different pathologic conditions.
- Sonora transducer hand pieces can be under water used.
- Wide range of preset protocols
- Customized programs
- Multi frequency hand pieces
- Automatic recognition of the connected accessory
- Continuous or pulsated delivery
- Associated with electrotherapy.

STANDARD FEATURES

- Sonora transducer hand pieces can be used underwater.
- Wide range of preset protocols
- Customized programs
- Multi frequency hand pieces
- Automatic recognition of the connected accessory
- Continuous or pulsated delivery
- Associated with electrotherapy

CURRENTS

- Diadynamic: MF MFS CP DF DFS CP
- Faradic: Rectangular Triangular Exponential Neo faradic Trabert
- Galvanic: Continuous MF Galvanic
- Impulsive: Rectangular Triangular Exponential Biphasic
- TENS: Burst Modulated Random
- Bipolar: Medium freq. Russian Bipolar IF Tetra polar Isoplanar Vectorial

SPECIFICATION

• Main input voltages	: 100-240V AC, 1.0 to 0.42A,
• Output frequency	: from 1 to 3 MHz
• Intensity	: 3 W/cm ²
• Output direct current)	: 30 mA
• Output pulse current	: 70 mA
• Output voltage	: 100 Vdc
• Power supply	: 90-240 Vac / 50-60 Hz

ORDERING DETAILS

CODE:	DESCRIPTION
SNR003	Sonora combination 2 Channel Electrotherapy & 1 & 3 Ultrasound