

IDEAL PLUS TWO SECTION COUCH

Designed and constructed using modern engineering techniques and manufactured from the finest material to the highest standard. Ideal Plus variable couches are the vision of tomorrow, featuring double gas struts configuration at head section, twin lift stability, open base design, self levelling feet and available in choice of 5 vibrant colours. The Ideal Plus range combines quality and durability with style and elegance and unrestricted access for a range of patient treatment and are supplied with breathing hole and paper towel holder as standard.




Painted / Chrome Side Rail




Single Side Padded Side Rail

Top Elevation:



Operational Parameters:



STANDARD FEATURES

- Twin lift design
- Positive/negative head section
- Retractable wheels and Paper Roll Holder as standard
- Breathing hole and plug as standard
- Electric Hi Lo operation by hand / foot control
- Sleek design with Rigid frame and rounded edges for safety
- Minimum height 46cm to maximum height 91cm
- Adjustable head section by dual gas struts from -20° to +80°
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm castors for safety and mobility
- Retractable dual locking castors on all 4 corners
- Padding thickness 5cm

SPECIFICATION

- Total Length : 195cm
- Width : 71cm
- Head section length : 68cm
- Leg section length : 124cm
- Height adjustable : 46cm to 91cm
- Loading weight : 225Kgs
- Upholstery : PVC vinyl with fire resistant
- Certification : CE

ORDERING DETAILS

CODE:	DESCRIPTION
TAB 200 / ELEC	2 SECTION ELECTRIC COUCH
TAB 200 / HYDR	2 SECTION HYDRAULIC COUCH
Optional Accessories : Padded / Painted SIDE RAIL & IV POLE HOLDER	

AVAILABLE COLORS



Manufactured By

