

IDEAL PLUS THREE SECTION COUCH



What definitely makes our Ideal Couch different is the pleasure in use and the true personality of the owners. Intensively tested for standard conformity, the couch meets the highest quality demands.

STANDARD FEATURES

- Positive/negative head and foot section
- Retractable wheels as standard
- Breathing hole and plug as standard
- Electric Hi Lo operation by hand / foot control
- Twin lift design with sleek frame and rounded edges for safety
- Minimum height 46cm to maximum height 91cm
- Adjustable head section by dual gas struts from -35° to + 85°
- Adjustable foot section by dual gas struts from -29° to +20°
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm castors for safety and mobility
- Padding thickness 5cm



Couch with Side Rail

SPECIFICATION

• Total Length	: 195 cm
• Width	: 65 cm
• Head section length	: 68 cm
• Middle section length	: 56 cm
• Leg section length	: 68 cm
• Height adjustable	: 46 cm to 91 cm
• Loading weight	: 225 kgs
• Upholstery	: PVC Vinyl with Fire Resistant
• Certification	: CE

ORDERING DETAILS

CODE:	DESCRIPTION
TAB 300 / ELEC	3 Section Electric Couch
TAB 300 / HYDR	3 Section Hydraulic