

IDEAL PLUS THREE SECTION BARIATRIC COUCH



An excellent solution for supporting Bariatric patients in both sitting and lying positions, a must for all Hospital Departments, Health Centres and Care Centres.

STANDARD FEATURES

- Positive/negative head and foot section
- Retractable wheels as standard
- Breathing hole and plug as standard
- Electric Hi Lo operation by hand / foot control
- Twin lift design with sleek frame and rounded edges for safety
- Minimum height 46cm to maximum height 91cm
- Adjustable head section by Electric motor from -35° to $+85^{\circ}$
- Adjustable foot section by Electric motor from -29° to $+20^{\circ}$
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm castors for safety and mobility
- Padding thickness 5cm

SPECIFICATION

• Total Length	: 195 cm
• Width	: 80 cm
• Head section length	: 68 cm
• Middle section length	: 56 cm
• Leg section length	: 68 cm
• Height adjustable	: 46 cm to 91 cm
• Loading weight	: 225 kgs
• Upholstery	: PVC Vinyl with Fire Resistant
• Certification	: CE

ORDERING DETAILS

CODE:	DESCRIPTION
TAB 300 / BAR / ELEC	3 Section Electric Couch
TAB 300 / BAR / HYDR	3 Section Hydraulic Couch

