

MAT TABLE

The Mat Therapy Tables are almost identical to the Bobath models in terms of specification and features but without the backrest facility. The single section top design provides a large and stable variable height support surface for a wide range of exercises and Bobath rehabilitation procedures.

SPECIFICATIONS

Safe working load and lifting capacity (from minimum height) of 250kg (550lbs=39 stone)
 Choice of upholstery width - 105cm (41") or 125cm (49") or extra large model 150cm (59")
 Choice of manual or electric height elevation
 Excellent height range from 43cm (17") to 97cm (37")
 Retractable wheel design and adjustable foot for stability on uneven surfaces
 Oilite® bearings on all pivot points
 Clearance beneath frame for mobile hoist
 Hand switch fitted as standard on electric models



All couches are covered by fire retardant material.
 Available colours are:



ORDERING DETAILS

CODE:

TAB 105
 TAB 125
 TAB 150

DESCRIPTION

MAT TABLE 105cm
 MAT TABLE 125cm
 MAT TABLE 150cm

IDEAL PLUS MULTI-PURPOSE TREATMENT

Ideal plus multi-purpose treatment chair with electric actuator height operation, gas strut backrest and detachable arm trays as standard.

ORDERING DETAILS

CODE:

TAB 600

DESCRIPTION

IDEAL PLUS MULTI-PURPOSE TREATMENT CHAIR

OPTIONAL EXTRAS

- SLA 25 Paper towel holder
- SLA 30 Optional foot switch electric models only
- SLA 35 Wrap around bar (side movement operation) for electric models only (powder coated finish)
- SLA 40 Wrap around bar (up / down movement operation) available in powder coated finish or chrome finish

"When placing orders please state upholstery colour and optional extras"

All couches are covered by fire retardant material.
 Available colours are:



SPECIFICATIONS

Minimum height	585mm/ 23.01"
Maximum height	1060mm/ 41.8"
Cushion width	640mm/ 25.2"
Head section length	740mm/ 29.14"
Mid section length	460mm/ 15.1"
Leg section length	480mm/ 18.9"
Total length	1740mm/ 58.5"
Head tilt angle	+90° to -20°
Mid tilt angle	+15° to 0°
Leg tilt angle	+30° to 70°
Lift mechanism	Electric
Breathing hole	1
Retractable wheels	Standard x 4
Weight rating	225kg / 496lb

