

IDEAL PLUS 6 SECTION DRIANAGE COUCH



The couch comes with many features as standard including a breathing hole and gas assisted head & foot rests. The postural drainage couch has a middle section assisted by gas springs to enable centre lift. This facilitates good drainage positions for all areas of the lung whilst comfortably supporting the patient's whole body.

SPECIFICATION

▪ Length	: 78"/198cm
▪ Width	: 27"/68.5cm 29"/73.6cm
▪ Electric Height Adjustment	: 46cm- 95cm (Dewert)
▪ Pneumatic Headrest Angle Adjustment	: 75°- 40°
▪ Pneumatic Leg rest Angle Adjustment	: 0°- 80°
▪ Foam Thickness	: 5cm
▪ Working Capacity	: 250kg
▪ Weight	: 79kg
▪ Accessories	: Removable Face Cushion, Regular Foot Controller

OPTIONS

- Head Section Design
- Front Arm rest
- Side Arm rest,
- Stand-by Castor System,
- Touch Stop Safety System,
- Paper Roll Holder

ORDERING DETAILS

CODE:	DESCRIPTION
TAB 601	6 Section Drianage Couch

IDEAL PLUS 7 SECTION ELECTRICAL COUCH



Ideal Plus 7 section couch incomparable armrests and gas strut head and foot sections. All other features are identical to the ideal plus 2 and 3 section couches.

SPECIFICATION

▪ Minimum height	: 560mm / 22.2"
▪ Maximum height	: 990mm / 38.9"
▪ Cushion width	: 700mm / 27.6"
▪ Head section length	: 390mm / 15.3"
▪ Leg section length	: 640mm / 33.0"
▪ Total length	: 1940mm / 76.3"
▪ Head tilt angle	: +25" to - 35"
▪ Lift mechanism	: Electric
▪ Breathing hole	: 1
▪ Retractable wheels	: 4
▪ Weight rating	: 225kg / 496lb

ORDERING DETAILS

CODE:	DESCRIPTION
TAB 700	7 Section Electrical Couch