

EIGHT - SECTION THERAPY COUCH



SPECIFICATION

- Table size : 65cm×184 cm
- Bed height range : 58-99 cm
- Neck section adjustable : Angle -60°- 40°
- Leg section adjustable : Angle 0-85°
- Middle section adjustable : Angle 0-130°
- User weight : 250 Kg
- Weight : 100Kg

FEATURES

- Eight Sections Design All. The Sections Can Be Adjusted
- Freely And Independently.
- Easy Access Foot Switch, 360° Accessible From Any Position Around The Table.
- Antibiosis With High Elasticity Sponge. Fire Resistance,
- Temperature Resistance, Durable And Comfortable.
- Stabilus Germany Air spring Which Makes The Operation Flexible And Safe.
- High Temperature Powder Coated Painting.
- The Middle Of The Table Is Controlled By Motor.

ORDERING DETAILS

CODE:	DESCRIPTION
TAB 800	Eight Section Therapy Couch

AVAILABLE COLORS

Black	Grey	Beige	Vivid Blue	Spruce Blue	Navy Blue

IDEAL PLUS 9 SECTION COUCH

FEATURES

- Nine sections and each section is adjustable independently.
- Adopt better motor with little working noise.
- The table is equipped with Germany "STABILUS" air spring which makes the operation flexible and safe.
- Equipped Collapsible and quiet truckle.
- 360° Foot control switch
- Antimicrobial and wear proof PU.
- Adjust range of handle standing 0-9mm
- Upright function hand remote control setting

SPECIFICATION

- Adjust angle of head : 0 to 9 tab to move
- Dimension : 2080 x 885 x 530-970mm
- Adjust angle of head : -20°-40°
- Adjust angle of back : 0°-25°
- Adjust angle of leg : 0-75°
- Loading bearing of head section : 127N
- Loading bearing of back section : 764N
- Loading bearing of leg section : 108N



ORDERING DETAILS

CODE:	DESCRIPTION
TAB 900-1	9 Section Electrical Couch