

IDEAL PLUS FOUR SECTION THERAPY COUCH



Ideal plus 4 section couch with head section and split leg section. Electric variable height table for drainage treatments, Designed to give maximum comfort and assistance in the positioning of the patient during the treatments. The provision of electric driven positioning permits effective and effortless positioning of even heavy patients for treatment of a number of pulmonary related conditions. Foot or Hand switch as standard.

SPECIFICATION

- Minimum height : 460mm / 18.1"
- Maximum height : 910mm / 35.8"
- Cushion width : 700mm / 27.6"
- Head section length : 360mm / 14.2"
- Mid section length : 510mm / 20.1"
- Leg section length : 920mm / 36.3"
- Total length : 1955mm / 77"
- Head tilt angles : 8 +30° to -50°
- Split leg tilt angles : +70° to -0.0°
- Mid section tilt : +13° to +20°
- Lift mechanism : Electric
- Breathing hole : 1
- Retractable wheels : Standard x 4
- Weight rating : 225kgs / 496lb

ORDERING DETAILS

CODE:	DESCRIPTION
TAB 401	4 Section Ideal Plus Couch

IDEAL PLUS FIVE SECTION COUCH



Designed and constructed using modern engineering techniques and manufactured from the finest material to the highest standard. Ideal plus variable couches are the vision of tomorrow, featuring double gas struts configuration at head section, twin lift stability, open base design, self leveling feet and available in choice of 5 vibrant colors, the Ideal Plus range combines quality and durability with style and elegance and unrestricted access for a range of patient treatment and are supplied with breathing hole and paper towel holder as standard.

SPECIFICATION

- Adjustable Height : 495mm / 19.5"
to 940mm / 37.0"
- Cushion width : 700mm / 27.6"
- Head section Length : 650mm / 25.6"
- Mid-Section Length : 680mm / 26.8"
- Leg Section Length : 660mm / 26"
- Total Length : 2060mm / 81.1"
- Head tilt angles : +87° to - 50°
- Leg Tilt Angles : +85° to - 50°
- Lift mechanism : Electric
- Breathing hole : 1
- Retractable wheels : Standard x 4
- Weight rating : 225kg / 496lb

ORDERING DETAILS

CODE:	DESCRIPTION
TAB 500	5 Section Electrical Couch

AVAILABLE COLORS



Black



Grey



Beige



Vivid Blue



Spruce Blue



Navy Blue