

ORTHOPEDIC ELECTRICAL EXAMINATION 7 SECTION COUCH



SPECIFICATION

- Height : Adjustable, Tilting, on Casters
- Operation : Electric
- No. of sections : 6-section
- Weight capacity : 200 kg
- Length : 185 cm
- Width : 69 cm
- Height: Min. : 60 cm
- Max. : 92 cm
- Mattress thickness : 5 cm

A variety of functions for most physiotherapy solutions. Quality assurance, good stability, silent operation, no disturbance, Excellent performance, evenly disperse the body's own weight, supporting all parts of the body, so that the muscles are completely relaxed. Tensile and bending strength has been rigorously tested. High-precision robot welding is enough to load static weight up to 175kg.

ORDERING DETAILS

CODE:	DESCRIPTION
TAB 700SL	7 Section Couch



IDEAL PLUS 7 SECTION ELECTRICAL COUCH



SPECIFICATION

- Length : 78"/198cm
- Width : 22"/56cm, 25"/64cm, 27"/69cm, 29"/74cm
- Height : 18"/46cm - 37"/95cm
- Foam : 2.5"/6.4cm
- Lifting Capacity : 440lbs/ 200kg
- Working Weight : 600lbs/272kg
- Product Weight : 174lbs/79kg

A very strong, robust and aesthetically pleasing treatment couch designed to meet the increasing need to cater for the heavy patient with electric height and back-rest functions as standard. Our seven section treatment couch features positioning casters with an adjustable self-leveling foot making them ideal for a wide range of healthcare environments and treatment procedures. Heavy duty rugged steel frame, durable and reliable actuator, Upholstered in fire proof PU/ leather, High density padding.

ORDERING DETAILS

CODE:	DESCRIPTION
TAB 710	7 Section Couch

