

THERMA-TECAR THERAPY

DIATHERMIC CAPACITIVE AND RESISTIVE TREATMENT



Therma is a temporary-use medical device, suitable for the application of therapeutic treatments diathermy capacitive and resistive energy transfer of musculoskeletal injuries (contractures, sprains and strains, tendonitis, sinusitis and bursitis) and osteo-articular pathologies, **Therma** is a new device for the treatment capacitive and resistive thermal noninvasive capable of generating, with the best efficacy and pain, hyperemia endogenous and stimulate the blood circulation in all the different cell layers to the lymphatic drainage. It has anti-inflammatory effects and results in an improvement of blood circulation and a significant reduction in pain. It stimulates a reaction of the natural and physiological body adequately to the specific pathology, and accelerates the healing process. Thanks to special applicators is able to provide the right amount of energy in the cell layers of interest. The therapeutic treatment thermal resistive and capacitive is the most recent and effective approach between diathermy treatments for all professionals who are looking for something more than a traditional diathermy unit (microwave or shortwave).



FEATURES

- Wide range of preset protocols
- Custom Programs
- Capacitive mode for drainage and de-contracting with greater thermal effect, resistive mode for biological stimulation with minimal thermal effect.
- Hand piece with 3 working configurations

SPECIFICATION

▪ Maximum power	: 150W / 250VA
▪ Working frequency	: 450 kHz / 500 MHz
▪ Operating mode	: Resistive / Capacitive
▪ Power	: 90-240 Vac / 50-60 Hz
▪ Max. power consumption	: 300 VA

ORDERING DETAILS

CODE:	DESCRIPTION
THR001	Therma Tecar Unit