

Intellect Transport 2 Combination Unit

The unit can be used with or without the optional cart, the lightweight, durable design and battery powered option provide the versatility needed to treat your patients and athletes in any professional environment. The Intellect Transport 2 is easy and intuitive to use and includes all of the essential features needed to get your patients moving.

All combo units feature Chattanooga's proprietary VMS waveform which enables stronger, more efficient muscle contractions along with 13 other types of electrical currents.

BETTER TREATMENT OUTCOMES

- Provides symptomatic relief of acute and chronic pain associated with musculoskeletal conditions and the management of post-operative pain1
- Improved ultrasound coupling detection and monitoring ensures the patient is receiving the intended energy for their treatment.
- Chattanooga's proprietary VMS waveform achieves stronger muscle contractions and less neural fatigue while still providing a more comfortable experience for the patient.
- Dual Frequency Ultrasound 1MHz or 3.3 MHz capability with 2cm2, 5cm2 and 10cm2 applicator options.

RELIABLE WITH IMPROVED USER EXPERIENCE

- Improved battery life removes the need for excess cables and provides all day use in portable applications.
- High resolution color screen with clear intuitive menus provide real time information on treatment progress.
- Device memory allows the end user to save up to 15 custom protocols for future use.
- Improved Responsiveness - Faster refresh rate and a more reactive key pad provides a better user experience for therapists and patients.
- Ultrasound Coupling Detection - The device is able to recognize whether or not the ultrasound applicator is making a good detection with the patient. When the treatment is interrupted the device will stop the treatment and/or beep to alert the therapist.

Technical Specification

Mains Power	100-240 VAC, 50/60 Hz, 1.0 A
Weight	2.2 kg (4.9 lbs)
Shipping Weight	3.9 kg (8.5 lbs)
Dimensions	34.6 cm x 25 cm x 12 cm (13.6"L x 9.8"W x 4.7"H)
Electrical Safety Class	Class 1, Electrotherapy Type BF, Ultrasound Type B
Safety Tests	IEC 60601-1, ISO 10993-1, IEC 60601-2-5, IEC 60601-1-2, IEC 60601-2-10, IEC 60601-1-6, IEC 62366-1, IEC 62304, IEC 62133, ISO 14971

Description	Part No.
Intellect Transport 2 Combination Unit	4778
System Cart	15-1136



WAVEFORMS

- VMS
- Interferential Current (Traditional IFC, 4-pole)
- Pre-modulated IFC (2-pole)
- High Voltage Pulsed Current
- Microcurrent
- Asymmetrical Biphasic current - TENS
- Symmetrical Biphasic current - TENS
- Galvanic
- Galvanic Interrupted
- Monophasic Triangular Pulsed
- Monophasic Rectangular Pulsed
- Träbert
- Russian
- Diadynamic LP, MF, DF, Cpid/Cpiso

Standard Accessories:

- 5cm2 Ultrasound treatment applicator,
- Electrode sponge
- 4 x carbon electrodes
- Stim ch 1/2 lead wire kit,
- 5cm round electrodes 10pk
- Ny LaTeX wrap (2 per pack)
- User manual , Power lead.