

LAMBDA YAG LASER 10W

Lambda Yag is a laser therapy equipment, which allows to work in depth on different types of pathologies. This, combined with the emission power, allows a good diffusion and the right absorption by human tissue, by combining an analgesic effect, anti-inflammatory effect and of the tissue regeneration. It allows to work in two modalities: continuous and pulsed. The continuous emission mode is suitable for pathologies where the injury is at depth. It is able to trigger cellular bio-stimulation processes, and the process of tissue regeneration, reducing movement recovery times. The pulsed emission mode is suitable for pathologies with acute phase, in order to perform anti-inflammatory pathologies with low thermal effects also suitable for the fast tissue regeneration and the maintaining of the cellular balance. This allows the controlled vasodilation and the better recovery of the inflammatory process.

ORDERING DETAILS

CODE:	DESCRIPTION
LMB003	Lambda Yag Laser up 10W With 7" Touchscreen



FEATURES

- Wavelength 1064nm
- Guide Light 650nm
- Laser power up to 10W
- Emission Continuous, Pulsed
- Power supply 90-240V/50-60Hz
- Absorption 100VA
- Laser class 4

INTELECT MOBILE LASER

- Class 3B Low Level Laser Therapy
- Independent control over all parameters
- Continuous and pulsed operation available
- Fully variable pulsed frequencies
- Real time feedback of energy delivered
- Display of energy in choice of joules or joules/cm
- Helpful menu of Clinical Indications specific to the laser applicator
- 15 user-defined memory slots on board
- Built in acupuncture point finder
- 17 different applicators available, including LED, laser diodes and clusters
- Innovative design allows for cart, table top, wall or mobile use
- Backlit LCD display for better visibility in low light settings
- Improved pq button for easier navigation.

LASER APPLICATORS (OPTIONALS)

• 27799 670nm	LED Diode 10mW
• 27803 850nm	Laser Diode 40mW
• 27840 850nm	Laser Diode 100mW
• 27804 850nm	Laser Diode 150mW
• 27841 850nm	Laser Diode 200mW
• 27805 820nm	Laser Diode 300mW



STANDARD ACCESSORIES

- Laser Protection Glasses (2x)
- Usermanual on CD

OPTIONAL ACCESSORIES

- Remote Interlock kit
- Intellect Mobile Carrying Bag
- Battery Pack
- Therapy System Cart
- Cart Adaptor

ORDERING DETAILS

CODE:	DESCRIPTION
CHA2779	Intellect Mobile Laser