

## YOUTECARX BASE READY FOR ACTION

YoutecarX can be fully customized by professionals for their own specific needs. There is also a choice of configurations designed for specific targets.



### STANDARD ACCESSORIES

- 1 power cable.
- 1 cable for unipolar applications.
- 1 unipolar hand piece.
- 1 steel neutral electrode.
- 1 resistive electrode Ø 61 mm.
- 1 capacitive electrode Ø 61 mm.
- Cable with clamp for disposable neutral electrode.
- Disposable gel adhesive electrode (5 pcs).
- Dimensions : 31cm x 16cm x 25cm
- Weight : 3.250 Kg

### CONFIGURATION

- Power : 150 Watt
- Frequency : 500KHz
- Output : Unipolar
- Emission : Continuous

### CORE & SYSTEM. SELECT POWER LEVEL, FREQUENCIES AND FORM OF EMISSION

#### CUSTOMIZATION OPTIONS

| CODE          | DESCRIPTION                                  |
|---------------|----------------------------------------------|
| • YTX-ELCAP41 | Capacitive Unipolar Electrode Ø 41 mm        |
| • YTX-ELCAP61 | Capacitive Unipolar Electrode Ø 61 mm        |
| • YTX-ELCAP81 | Capacitive Unipolar Electrode Ø 81 mm        |
| • YTX-ELRES41 | Resistive Unipolar Electrode Ø 41 mm         |
| • YTX-ELRES61 | Resistive Unipolar Electrode Ø 61 mm         |
| • YTX-ELRES81 | Resistive Unipolar Electrode Ø 81 mm         |
| • ELNH        | Steel Neutral Plate 14x20 cm                 |
| • CAVO-U      | Cable for Unipolar Applications Res/Cap      |
| • YTX-MAN-U   | Hand piece for Unipolar Applications Res/Cap |

### ORDERING DETAILS

| CODE:  | DESCRIPTION          |
|--------|----------------------|
| YTX001 | YoutexCarx Base Unit |



**POWER**  
150 Watt (included)  
200 Watt  
250 Watt



**FREQUENCY**  
500 KHz (included)  
Standard Pack 3 freq.  
(500 KHz / 750 KHz / 1 MHz)  
On Demand Pack



**EMISSION**  
Continuous (included)  
Sensitive Athermy  
HEA (High Energy Athermy)  
Athermy Pack (Sensitive + HEA)

**BIPOLAR** - Work without neutral plate using BIPOLAR accessories Requires implementation of the AUX outlet.

| CODE         | DESCRIPTION                                                           |
|--------------|-----------------------------------------------------------------------|
| YTX-ELCAP41B | Capacitive Bipolar Electrode Ø 41 mm                                  |
| YTX-ELCAP61B | Capacitive Bipolar Electrode Ø 61 mm                                  |
| YTX-ELCAP81B | Capacitive Bipolar Electrode Ø 81 mm                                  |
| YTX-VERT     | Resistive Bipolar Paravertebral Electrode (adjustable heads distance) |
| YTX-MAN-B    | Hand piece for Res/Cap Bipolar Applications                           |



**NEUTRODINAMIC** Dynamic applications with two active hand pieces and without neutral plate, using NeuroDinamico® accessories. Requires implementation of the AUX outletO®.

| CODE       | DESCRIPTION                      |
|------------|----------------------------------|
| CAVO-ND    | Multifunctional Aux Cable        |
| YTX-MAN-ND | NeuroDinamico® Hand piece        |
| YTX-ELND41 | NeuroDinamico® Electrode Ø 41 mm |
| YTX-ELND61 | NeuroDinamico® Electrode Ø 61 mm |
| YTX-ELND81 | NeuroDinamico® Electrode Ø 81 mm |



**AUTOMATIC APPS** - Not operator-dependent Requires implementation of the AUX outlet.

**ELAUTO** Kit for Static Automatic Applications, including: Main cable and connections, 4 resistive clips electrodes. 2 elastic bands.

**JOULE MODE** Set your treatment by controlling the energy level in Joules. Requires implementation of the AUX outlet.

#### DESCRIPTION

- Enables the Joule Mode Software
- Includes free Contour Automatic Applications Kit:
- 8 Contour electrodes large size
- 8 Contour electrodes small size
- 2 Contour electrodes neutral size
- 1 Black Hub
- 1 White Hub
- 2 Elastic Bands 100 cm



**INTIMITY** - Pelvic Floor rehabilitation Requires implementation of the AUX outlet.

| CODE         | DESCRIPTION                                                                                                                |
|--------------|----------------------------------------------------------------------------------------------------------------------------|
| ENDO         | Endoprobe Unipolar vaginal Probe Ø 20 mm                                                                                   |
| ENDOB        | Endoprobe Bipolar vaginal Probe Ø 20 mm                                                                                    |
| ENDOA        | Endoprobe Unipolar anal Probe Ø 12 mm                                                                                      |
| YTX-MAN-B    | Handpiece for Res/Cap Bipolar Applications                                                                                 |
| YTX-ELCAP41B | Capacitive Bipolar Electrode Ø 41 mm                                                                                       |
| ELAUTO       | Kit for Static Automatic Applications including: Main cable and connections, 4 resistive clips electrodes, 2 elastic bands |
| POINT        | Unipointer                                                                                                                 |



#### CREAM & GEL:

Youtecar Cream and Gel rely on brand new formulas, specifically prepared to have the maximum effectiveness



#### TROLLEY:

The YTX trolley is specifically designed to contain the device and all the accessories! (IATA Cabin Size)



#### CART:

With the YoutecarX cart you can keep all your accessories on hand

